

AU5800 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Reagents								
• Check reagent status	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Monitor the Reagent Status ↳ Monitor the Reagent Status ↳ Replace the Reagents	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Load reagent bottles in the refrigerators		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform a reagent check		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Check ISE reagent status	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Perform the ISE Startup ↳ Inspect the ISE Reagents ↳ Replace the ISE Reagents	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Load ISE reagents		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Assign a reagent to fixed position	AU5800 Instructions for Use Chapter 5 System Monitoring and Results ↳ Reagent Management ↳ Assign a Reagent Position	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Identify reagent preparation	Reagent Instructions for Use ↳ Reagent Preparation	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Calibration								
• Determine analyzer calibration status	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Monitor the Reagent Status	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform analyzer calibration	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Calibrate Tests	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform ISE calibration	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Perform the ISE Startup ↳ Calibrate the ISE	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

AU5800 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Process QC								
• Perform daily QC	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Process Quality Control (QC)	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Order QC for multiple bottles (Individual Requisition)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Order and process QC for specific tests		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Process Patient Samples								
• Identify correct sample racks and containers for sample processing	AU5800 Instructions for Use Chapter 4 Sample Programming and Processing ↳ Sample Preparation ↳ Place the Sample Cups or Tubes in the Rack ↳ Prepare Racks for Analysis ↳ Adding Racks Directly to the Trays on the Rack Input Component ↳ Adding Racks to the Priority Rack Input Component ↳ Order (Requisition) for Routine and Emergency Samples	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Process bar coded samples from the rack input component		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Process samples using the priority rack input component		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Manually order samples for processing on routine racks		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Manually order samples for processing on emergency racks		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Review Results								
• Review Calibration Monitor	AU5800 Reference Manual Chapter 4 System Monitoring and Results ↳ Monitor the Reagent Blank and Calibration ↳ Monitor QC	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Review QC Monitor		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

AU5800 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Review Results, continued								
• Monitor Sample Status	AU5800 Instructions for Use Chapter 5 System Monitoring and Results ↳ Monitoring Analysis ↳ Sample Status Screen	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Print results	AU5800 Instructions for Use Chapter 4 Sample Programming and Processing ↳ Print Results	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Identify flags	AU5800 Instructions for Use Chapter 7 Flags	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Maintenance								
• Perform Daily Maintenance (ISE and Analyzer)	AU5800 Instructions for Use Chapter 6 Maintenance ↳ Daily Maintenance ↳ Weekly Maintenance ↳ Monthly Maintenance ↳ Quarterly Maintenance ↳ Six-Month Maintenance ↳ ISE Maintenance for All Markets Except Japan	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform Weekly Maintenance (ISE and Analyzer)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform ISE Every Other Week Maintenance (ISE only)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform Monthly Maintenance (Analyzer only)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform ISE Every Other Month Maintenance (ISE only)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform Quarterly Maintenance (ISE and Analyzer)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform 6-Month Maintenance (ISE and Analyzer)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

AU5800 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
Maintenance, continued								
<ul style="list-style-type: none"> Perform Yearly Maintenance (Analyzer only) 	AU5800 Instructions for Use Chapter 6 Maintenance ↳ Yearly Maintenance	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
Additional Procedures								
<ul style="list-style-type: none"> Set a new index 	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Set a New Index	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> Check Analyzer Status 	AU5800 Instructions for Use Chapter 2 Daily Startup ↳ Inspect the Analyzer Status	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> Perform an End Process 	AU5800 Instructions for Use Chapter 5 System Monitoring and Results ↳ System Shutdown (End Process)	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> Perform and recover from an Emergency Stop 	AU5800 Instructions for Use Chapter 5 System Monitoring and Results ↳ Perform an Emergency Stop	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> 		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> 		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> 		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> 		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			