

AU5800 Monthly Maintenance Job Aids



For Training Purposes Only

These job aids are shortened versions of the procedures found in the source below. The procedures are written as standalone procedures to ensure they can be performed in any order. The optional ISE module has references for two flow cells. Information in the job aid is correct as of the date published. Verify you have the correct information. Part numbers for supplies required for the maintenance procedures are found in the reference listed below.

Source: AU5800[®] Chemistry Analyzer User's Guide PN A98352AB (October 2012)


Document Disclaimers

Document Disclaimers	This document is not intended to replace the information in your User's Guides, Quick Response Guide or other product documentation. Information in the User's Guide and Quick Response Guide supersedes information in any other manual.
Warnings and Precautions	<p>Read all product manuals and consult with Beckman Coulter trained personnel before attempting to operate the instrument.</p> <p>Beckman Coulter, Inc. urges its customers and employees to comply with all national health and safety standards such as the use of barrier protection. This may include but is not limited to, protective eyewear, gloves, suitable laboratory attire when operating or maintaining this or any other automated laboratory equipment.</p>
Intention for Use	This document is not intended to replace the information in your User's Guide and Quick Response Guide. Information in the User's Guide supersedes information in any other manual.
Revision Status	Rev. A (May 2013) Software version 4.11
Trademarks	AU5800® Chemistry Analyzer


Clean the Sample Probe and Reagent Probe Wash Wells

Supplies Required:

- Cotton swabs
- Disposable transfer pipette
- 10% Bleach (Beckman Sodium Hypochlorite (5%) diluted 1:10 concentration. Prepare by adding 10 parts bleach to 90 parts DI water)



Clean the Sample Probe and Reagent Probe Wash Wells continued page 2



Clean the Mix Bar Wash Wells

- Supplies Required:**
- Cotton swabs
 - Disposable transfer pipette
 - 10% Bleach (Beckman Sodium Hypochlorite (5%) diluted 1:10. Prepare by adding 10 parts bleach to 90 parts DI water)

Ensure the system is in **Standby Mode**

Open the front and rear main cover on each analyzer unit

Manually turn the mix bar units so the mix bars are away from the wash wells

Use a transfer pipette to dispense 10% bleach solution into each wash well. Use a separate cotton swab to clean each wash well. Avoid splashing and clean up spills immediately

From the Home screen, select the **Analyzer Maintenance** jump button

Place a ✓ in the **Analyzer Maintenance** check box

Select **Replace Mix Bar**. Make the following selections:

- All units are selected by default, deselect units not required
- Select **The First Mixer**
- At **Times**, enter 1, select **OK**

Press the green **DIAG** button to initialize the R1/S mix bar unit

Select **Replace Mix Bar**. Make the following selections:

- All units are selected by default, deselect units not required
- Select **The Second Mixer**
- At **Times**, enter 1, select **OK**

Press the green **DIAG** button to initialize the R2 mix bar unit

Select **Prime Washing Line**. Make the following selections:

- All units are selected by default, deselect units not required
- At **Times**, enter 1, select **OK**

Press the green **DIAG** button. The analyzer will prime DI water through the sample and reagent probes, mix bars, and the wash nozzle. Inspect the wash wells for proper drainage


Are wash wells draining properly?

Document you completed the procedure on the paper maintenance log

Deselect the **Analyzer Maintenance** check box. Close the front and rear main covers

Clean the Wash Nozzle Unit and Check the Tube Mounting Joints

- Supplies Required:**
- Dry clean cloth or paper towel
 - Sonicator with DI water



Clean the Wash Nozzle Unit and Check the Tube Mounting Joints continued page 2

Inspect and clean each wash nozzle unit by placing it in a sonicator filled with DI water for 15 minutes. Only the nozzles need to be submerged. Do not wet or place the springs in the DI water. If the springs get wet, wipe them with a dry clean lint free cloth

Remove each wash nozzle unit from the sonicator and dry thoroughly with a dry clean lint free cloth

Install the wash nozzle unit by placing it over the positioning screws and tighten the silver knob to hold it in place

Inspect the packing (thin ring) for each manifold to ensure it is not worn or damaged. Install each manifold using the color coded positions for proper placement. Tighten the manifolds finger tight to prevent cuvette overflow but do not over-tighten

Yes

Do you have another wash nozzle unit to install?

No

Select **Prime Wash Nozzle**, (systems with multiple units are all selected by default, deselect units if not required), for **Times** enter 5, select **OK**

Press the green **DIAG** button. Verify the wash nozzle unit moves up and down without interference and that there are no leaks

Deselect the **Analyzer Maintenance** check box

Close the rear main cover of each analyzer unit

Document you completed the procedure on the paper maintenance log

Note: A sonicator is recommended for cleaning the nozzles but if one is not available, clean the nozzles with the supplied stylet and DI water

Clean the DI Water Tank, DI Water and Sample Probe Filters

- Supplies Required:**
- Dry clean lint free cloth
 - Basin
 - 20% Bleach (Beckman Sodium Hypochlorite (5%) diluted 1:20. Prepare by adding 20 parts bleach to 80 parts DI water)
 - DI water tank filled with 5 L of DI water
 - Sonicator with DI water

Perform this procedure for each analyzer unit. Ensure the system is in **Standby**. From the Home screen, select the **End Process** button

Open the front right door on each analyzer unit

Place a basin on the floor in front of the DI Water tank to catch spilled water

Disconnect the float sensor connector #795

Press the grey quick disconnect joints on the front of the tank, remove the tubings from the front of the tank

Pull the DI water tank forward and pull out the tubings from the top of the tank. Wrap the tubings with a dry clean lint-free cloth (1 tube has the DI water filter attached)

Pull the DI water tank out of the analyzer. Unscrew and remove the float sensor from the DI water tank, wipe the float sensor with DI water, then wipe the float sensor dry with a dry clean lint free cloth. Install the float sensor in the spare DI water tank containing 5 L of DI water

Discard the DI water from the tank you removed from the analyzer and rinse the tank with the 20% bleach solution. Rinse the tank thoroughly with DI water and allow to dry

Locate the filter on the DI water supply tubing removed from the tank. Position the tubing over the basin and unscrew the filter case, remove the filter and allow the water to drip

Be sure to remove the float sensor connector to prevent water from pouring out of the DI water tubing

Avoid pulling the DI water tank by the front valves

Continue on the next page

Clean the DI Water Tank, DI Water and Sample Probe Filters continued page 2

