

## AU680 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
<b>Reagents</b>								
• Check reagent status	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Monitor the Reagent Status ↳ Replace the Reagents	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Load reagent bottles in the refrigerators		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform a reagent check		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Check ISE reagent status	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Perform the ISE Startup ↳ Inspect the ISE Reagents ↳ Replace the ISE Reagents	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Load ISE reagents		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Assign a reagent to fixed position	<a href="#">AU680 Instructions for Use</a> Chapter 5 System Monitoring and Results ↳ Reagent Management ↳ Assign a Reagent Position	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Identify reagent preparation	<a href="#">Reagent Instructions for Use</a> ↳ Reagent Preparation	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<b>Calibration</b>								
• Determine analyzer calibration status	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Monitor the Reagent Status	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform analyzer calibration	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Calibrate Tests	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform ISE calibration	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Perform the ISE Startup ↳ Calibrate the ISE	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

## AU680 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
<b>Process QC</b>								
• Perform daily QC	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Process Quality Control (QC)	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Order QC for multiple bottles (Individual Requisition)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Order and process QC for specific tests		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<b>Process Patient Samples</b>								
• Identify correct sample racks and containers for sample processing	<a href="#">AU680 Instructions for Use</a> Chapter 4 Sample Programming and Processing ↳ Sample Preparation ↳ Place the Sample Cups or Tubes in the Rack ↳ Placing Racks on the Rack Supply Component ↳ Order (Requisition) for Routine and Emergency Samples	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Process bar coded samples in racks		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Manually order samples for processing in racks		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Process bar coded samples in the STAT table	<a href="#">AU680 Instructions for Use</a> Chapter 4 Sample Programming and Processing ↳ Priority STAT Samples	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Manually order sample for processing in the STAT table		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<b>Review Results</b>								
• Review Calibration Monitor	<a href="#">AU680 Reference Manual</a> Chapter 4 System Monitoring and Results ↳ Monitor the Reagent Blank and Calibration ↳ Monitor QC	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Review QC Monitor		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

## AU680 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
<b>Review Results, continued</b>								
• Monitor Sample Status	<a href="#">AU680 Instructions for Use</a> Chapter 5 System Monitoring and Results ↳ Monitoring Analysis ↳ Sample Status Screen	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Print results	<a href="#">AU680 Instructions for Use</a> Chapter 4 Sample Programming and Processing ↳ Print Results	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Identify flags	<a href="#">AU680 Instructions for Use</a> Chapter 7 Flags	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<b>Maintenance</b>								
• Perform Daily Maintenance (ISE and Analyzer)	<a href="#">AU680 Instructions for Use</a> Chapter 6 Maintenance ↳ Daily Maintenance ↳ Weekly Maintenance ↳ Monthly Maintenance ↳ Quarterly Maintenance ↳ Six-Month Maintenance ↳ ISE Maintenance for All Markets Except Japan	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform Weekly Maintenance (ISE and Analyzer)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform ISE Every Other Week Maintenance (ISE only)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform Monthly Maintenance (Analyzer only)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform ISE Every Other Month Maintenance (ISE only)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform Quarterly Maintenance (ISE and Analyzer)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
• Perform 6-Month Maintenance (ISE and Analyzer)		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			

## AU680 Chemistry Analyzer

Task	Resources	How Measured?	Achieved?		N/A	Date	Trainee Initial	Trainer Initial
			Y	N				
<b>Maintenance, continued</b>								
<ul style="list-style-type: none"> <li>Perform Yearly Maintenance (Analyzer only)</li> </ul>	<a href="#">AU680 Instructions for Use</a> Chapter 6 Maintenance ↳ Yearly Maintenance	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<b>Additional Procedures</b>								
<ul style="list-style-type: none"> <li>Set a new index</li> </ul>	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Set a New Index	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li>Check Analyzer Status</li> </ul>	<a href="#">AU680 Instructions for Use</a> Chapter 2 Daily Startup ↳ Inspect the Analyzer Status	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li>Perform an End Process</li> </ul>	<a href="#">AU680 Instructions for Use</a> Chapter 5 System Monitoring and Results ↳ System Shutdown (End Process)	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li>Perform and recover from an Emergency Stop</li> </ul>	<a href="#">AU680 Instructions for Use</a> Chapter 5 System Monitoring and Results ↳ Perform an Emergency Stop	<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li></li> </ul>		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li></li> </ul>		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li></li> </ul>		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li></li> </ul>		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			
<ul style="list-style-type: none"> <li></li> </ul>		<input type="checkbox"/> Direct Observation <input type="checkbox"/> Test <input type="checkbox"/> Other _____	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			