


DxC 700 AU Daily Startup Using the STAT Table with Automation



DxC 700 AU Daily Startup Using the STAT Table with Automation, continued

✓ Perform the ISE Startup (for labs with an ISE Module)

1. Inspect the ISE reagents, and replace reagents that will expire onboard remaining or have a volume below daily requirements
 - Select **REAGENT > Reagent Management > ISE**
 - To replace a reagent:
 - Select the ISE reagent to replace > select **Replace ISE Reagent (F1)**
 - Confirm the ISE operation LED is off
 - Scan the bar code label on the bottle with the hand scanner or type in the lot number and expiration date
 - Replace the reagent bottle and select **OK**
 - Select **ISE Maintenance >** select the **ISE Maintenance** box
 - If ISE Buffer Solution is replaced, select **Buffer Prime >** select **OK**
 - If ISE MID Standard or ISE Reference Solution is replaced, select **MID/REF Prime >** select **OK**
 - Press the **green DIAG/TABLE ROTATION** button to start the prime
 - Clear the **ISE Maintenance** box when the prime is complete
2. Clean the ISE
 - Select **MAINT. > ISE Maintenance > Maintenance**
 - Fill a Hitachi cup with 1 mL of ISE Cleaning Solution and place the cup in the **CLEAN** position on the STAT table
 - Select **Cleaning**, select **OK** to begin the clean
 - When cleaning is complete, remove and discard the Hitachi cup
3. Perform a Total Prime if calibrating immediately following an ISE Clean
 - Select **MAINT. > ISE Maintenance > Maintenance >** select the **ISE Maintenance** box
 - Select **Total Prime**, select **OK**
 - Press the **green DIAG/TABLE ROTATION** button to start the prime
 - Clear the **ISE Maintenance** box when the prime is complete
4. Calibrate the ISE
 - Place Hitachi cups filled with approximately 500 µL of the required ISE Serum and/or Urine Standards into the **S-L, S-H, U-L, and U-H** positions of the STAT table
 - Select **Calibration** tab from ISE Maintenance
 - Select **Serum Start, Urine Start, or Serum/Urine Start**, select **OK** to begin the ISE calibration
 - When the calibration is complete, confirm the results are within the Slope and MID Factor ranges for serum and/or urine
 - Remove and discard the Hitachi cups
5. Remove the Operator Name: Select **HOME > Change Operator (F2)**, remove name or initials, select **OK**

✓ Order Calibrations and QC

1. Select **STAT > STAT (Calibration)**
 - In **Type**, select the sample type
 - Select **Auto CAL/QC Order (F4)**. From the STAT (Calibration) dialog, select **OK**. The system displays the reagent blank and calibrators ordered in blue. Confirm the ordered tests are correct
 - Select **Type**, to verify orders of all sample types in use
2. Select the **STAT (QC)** button
 - In **Type**, select the sample type
 - Select **Select All Tests (F5)**, select individual tests or select **Panel** to select a panel, and then select **OK**
 - Repeat for each sample type in use
5. Select **Save (F2)**
6. Select **CAL/QC Position (F7)**

✓ Prepare Reagent Blank, Calibrators and QC

Use the CAL/QC Position list to prepare the required reagent blank, calibrator and/or QC samples.
Use the list to ensure you use the correct bar coded tubes and required volume
(dead volume of the specific container is not included) for each sample. Select **Close**


Continue on the next page

DxC 700 AU Daily Startup Using the STAT Table with Automation, continued

Continued from the previous page

✓ Load RB/CAL/QC on the STAT Table

Note: The reagent blank, calibrators and QC can be loaded on the STAT table by opening the small STAT table cover or removing the large STAT table cover. If there are more than 22 calibrator and/or QC samples required in the run, load the additional calibrator and QC samples by removing and replacing the calibrators and QC using the same STAT table or by removing and replacing the entire STAT table with another prepared STAT table.


- Load the calibrators and/or control samples (maximum of 22) with the bar codes facing out into the STAT table. Calibrators must be loaded before QC. The reagent blank must be loaded in the **RB1** position
- Select  from the action button area

- From the STAT Start dialog, review errors on the Error List and perform any corrective actions, if necessary
- Select the **Last Table** box in the STAT Start dialog
- Select **Start**

When the STAT TABLE ROTATION LED stops blinking the aspiration is complete. There will be an audible event displayed stating aspiration is complete. Remove all the calibrator and/or QC samples

Is this the last load of calibrators and/or QC samples?

- From the STAT Start dialog, review errors on the Error List and perform any corrective actions, if necessary
- Select **Start**

- When the STAT TABLE ROTATION LED stops blinking, add more calibrators and/or QC samples
- Select  from the action button area

✓ Review RB/CAL/QC Results

1. Review the reagent blank, calibration, and QC reports for flags. Take appropriate actions based on flags
2. Verify reagent blank and calibration meet lab requirements. **MENU > Calibration > Calibration Monitor > Status**, select **Reagent Blank** or **Calibration** column for the desired test
3. Verify QC results meet lab requirements. **QC > Main** tab (all tests are selected by default) > **Chart View**. Review QC test results as required

Note:

Take appropriate steps at your connected automation computer to **Resume** sample loading to the DxC 700 AU.