# HEMATOLOGY PERFORMANCE VERIFICATION MANUAL

4277076AC July 2018

## TRADEMARKS

Beckman Coulter, the stylized logo, and the Beckman Coulter product and service marks mentioned herein are trademarks or registered trademarks of Beckman Coulter, Inc. in the United States and other countries.

All other trademarks, service marks, products, or services are trademarks or registered trademarks of their respective holders.

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFICATION MANUAL (mm/dd/yy) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

# **REVISION HISTORY**

This document applies to integrating your new Hematology System into your laboratory. When subsequent verifications affect the information in this document, a new issue will be released to the Beckman Coulter Web site. For labeling updates, go to www.beckmancoulter.com and download the latest version of the manual for your instrument.

Revision AA Initial Issue (August 2009)

Revision AB (July 2011)

Revision AC (July 2018)

The following sections were modified:

- Revision update of document from AB (July 2011) to AC (July 2018).
- Formatted procedures throughout the document to SOP guidelines.
- Page v of x: Safety Notice © 2011 Beckman Coulter, Inc. All Rights Reserved changed to © 2018 Beckman Coulter, Inc. All Rights Reserved.
- Page vii of x: Hematology Customer Support Contacts updated to current contact information.
- Page ix of x: Table of Contents updated to reflect chapters within the manual.
- Page 1-3 of 10: Deleted Prerequisites from the chart.
- Page 1-5 of 10: Trained Operator Responsibilities, first bullet changed word from reflex criteria to decision rule criteria.
- Page 1-6 of 10: Performed for Implementation and periodically (as per lab protocol and/or local regulatory agency), first bullet – added the word verify. Tenth bullet – added word Analytical and AMR.
- Assistance from your Applications Specialist, forth bullet added word Analytical and AMR.
- Page 1-7 of 10: Updated Reference at bottom of page.
- Page 2-9 of 10: Added section Reproducibility/Repeatability.
- Page 6-2 of 14: Important Note added System Messages.
- Page 7-17 of 28: Establishing Lab Limits Considerations, updated the web address <u>http://www.jointcommission.org</u> and date.
- Page 7-18 of 28: Updated CAP Requirements: Hematology from Coagulation Checklist: 09/27/2007 to Coagulation Checklist: 07/28/2015.
- Page 7-20 of 28: On the Note, corrected A@ to A2. Updated the References at the bottom of the page.
- Page 7-23 of 28: Automated Mixing Systems, added DxH Series to the first paragraph. Reworded the last paragraph in the section.
- Page 7-25 of 28: Updated References at the bottom of the page.
- Page 10-3 of 4: Clinical Laboratory Standards Institute (CLSI) Guidelines, updated the Standards and Guidelines.
- Deleted Chapter 11 Graphs.

# SAFETY NOTICE

READ ALL PRODUCT MANUALS AND CONSULT WITH BECKMAN COULTER-TRAINED PERSONNEL BEFORE ATTEMPTING TO OPERATE INSTRUMENT. DO NOT ATTEMPT TO PERFORM ANY PROCEDURE BEFORE CAREFULLY READING ALL INSTRUCTIONS. ALWAYS FOLLOW PRODUCT LABELING AND MANUFACTURER'S RECOMMENDATIONS. IF IN DOUBT AS TO HOW TO PROCEED IN ANY SITUATION, CONTACT YOUR BECKMAN COULTER REPRESENTATIVE.

#### HAZARDS AND OPERATIONAL PRECAUTIONS AND LIMITATIONS

WARNINGS, CAUTIONS and IMPORTANTS alert you as follows:

WARNING	-	Can cause injury.
CAUTION	-	Can cause damage to the instrument.
IMPORTANT	-	Can cause misleading results.

BECKMAN COULTER, INC. URGES ITS CUSTOMERS TO COMPLY WITH ALL NATIONAL HEALTH AND SAFETY STANDARDS SUCH AS THE USE OF BARRIER PROTECTION. THIS MAY INCLUDE, BUT IT IS NOT LIMITED TO, PROTECTIVE EYEWEAR, GLOVES AND SUITABLE LABORATORY ATTIRE WHEN OPERATING OR MAINTAINING THIS OR ANY OTHER AUTOMATED LABORATORY ANALYZER.

#### WARNING

Risk of operator injury if:

- All doors, covers and panels are not closed and secured in place prior to and during instrument operation.
- The integrity of safety interlocks and sensors is compromised.
- Instrument alarms and error messages are not acknowledged and acted upon.
- You contact moving parts.
- You mishandle broken parts.
- Doors, covers and panels are not opened, closed, removed and/or replaced with care.
- Improper tools are used for troubleshooting.

To avoid injury:

- Keep doors, covers and panels closed and secured in place while the instrument is in use.
- Take full advantage of the safety features of the instrument. Do not defeat safety interlocks and sensors.
- Acknowledge and act upon instrument alarms and error messages.
- Keep away from moving parts.
- Report any broken parts to your Beckman Coulter Representative.
- Open/remove and close/replace doors, covers and panels with care.
- Use the proper tools for troubleshooting.

#### CAUTION

System integrity might be compromised and operational failures might occur if:

- This equipment is used in a manner other than specified. Operate the instrument as instructed in the Product Manuals.
  - You introduce software that is not authorized by Beckman Coulter into your computer. Only operate your system's computer with software authorized by Beckman Coulter.
  - You install software that is not an original copyrighted version. Only use software that is an original copyrighted version to prevent virus contamination.

#### IMPORTANT

If you purchased this product from anyone other than Beckman Coulter or an authorized Beckman Coulter distributor, and, if it is not presently under a Beckman Coulter service maintenance agreement, Beckman Coulter cannot guarantee that the product is fitted with the most current mandatory engineering revisions or that you will receive the most current information bulletins concerning the product. If you purchased this product from a third party and would like further information concerning this topic, call your Beckman Coulter Representative.

© 2018 Beckman Coulter, Inc. All Rights Reserved

# **HEMATOLOGY CUSTOMER SUPPORT CONTACTS**

SALES AND SERVICE OFFICES FOR FIELD SERVICE OR TECHNICAL INFORMATION

CUSTOMER TECHNICAL SUPPORT C United States Canada	ENTER / SERVICE 800 526-7694 800 526-7694
ORDER ASSISTANCE United States	
BECKMAN COULTER Web Site	www.BeckmanCoulter.com
Sales Office: 250 South Kraemer Boulevard Brea, CA. 92821 800 526-3821	
Instrument Model	
Serial Number	
SID#/Instance#	
IQAP Participant #	
Your Account #	
Contacts:	
Customer Service Representative (Order/Contract Assistance)	
Applications Specialist	
Field Service Engineer	
Instrument Sales Representative	

# TABLE of CONTENTS

**Revision History, iii** 

Safety Notice, v

Hematology Customer Support Contacts, vii

INTRODUCTION, 1-1 Beckman Coulter Hematology Analyzer, 1-1 Performance Verification Flow chart, 1-3 Performance Verification Studies Overview, 1-4 Verification Synopsis, 1-5 Performance Verification Checklist, 1-9

PRE-CALIBRATION, CALIBRATION & QUALITY ASSESSMENT, 2-1

**CBC PARAMETER COMPARISON, 3-1** 

**DIFF PARAMETER COMPARISON, 4-1** 

**RETIC PARAMETER METHOD COMPARISON, 5-1** 

**TRUTH TABLE ANALYSIS, 6-1** 

OTHER STUDIES, 7-1

BODY FLUIDS COMPARISON, 7-3 MEASURING RANGE (LINEARITY), 7-7 MODE TO MODE COMPARISON, 7-9 ESTABLISH LAB LIMITS FOR QC, 7-15 REF INTERVAL STUDY (Normal Range), 7-21 MIXING STUDIES, 7-23

COMMENTARY, 8-1 Manual Diff versus Automated Diff, 8-1

GLOSSARY, 9-1

**REFERENCES**, 10-1

1

#### INTRODUCTION

## **Beckman Coulter Hematology Analyzer**

The Hematology Verification Manual is designed to assist you in making a smooth transition integrating your new Hematology System into your laboratory. Please take the time to review this book thoroughly.

Within this Performance Verification Manual you will find guidelines for evaluating your new analyzer. You will find procedures and worksheets to assist you with each stage of the verification process.

Please note that the procedures provided outline the general steps for characterizing the performance of your new system. Your current laboratory policies or your local regulatory agency may dictate more specific procedures for your laboratory to follow.

We suggest that you use this Hematology Performance Verification Manual to store all of your installation records and subsequent characterization data for you system.

We hope you will find the information contained in this manual a useful tool for making a smooth transition to your new hematology system.

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFICACTION (MANUAL(mm/dd/yy)) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released



# PERFORMANCE VERIFICATION STUDIES OVERVIEW

### CLIA '88

Much has been said and written about the impact of the Clinical Laboratory Improvement Amendments of 1988 (CLIA'88) on the laboratory and method evaluation.<sup>1</sup> It is often confusing to determine what is needed for accreditation and certification for the different types of testing categories.

The Centers for Medicare & Medicaid Services (CMS) regulates all laboratory testing (except research) performed on humans in the U.S. through the Clinical Laboratory Improvement Amendments (CLIA).

The following procedures will assist you in the evaluation of a new instrument or clinical test. Not all of these procedures are required by CLIA; however, some of them may be required by the College of American Pathologists (CAP), the Joint Commission of Accreditation of Healthcare Organizations (known as JCAHO), and/or individual agencies. It is important to know local requirements as well as those of any other accrediting agencies that will impact your laboratory. The data generated from some of these procedures will help to establish the baseline characterization of your system.

Your laboratory must compile its own policies and procedures manual for method evaluation, in compliance with the appropriate accrediting agencies. You are empowered to make your own decisions as to what procedures are appropriate, and which performance limits or specifications are acceptable.

# VERIFICATION SYNOPSIS

### Hardware Installation

The instrument installation will be performed by a local Beckman Coulter service representative, who will perform various system checks to ensure and document that the system meets specific Beckman Coulter performance specifications. This process should take between one and two days. You will be contacted to schedule this installation.

### **Instrument Implementation**

For the applications and performance verification process of your new analyzer, a Beckman Coulter Applications Specialist will assist you in using this Performance Verification Manual. An Applications Specialist will contact you to schedule time with the assigned Trained Operator.

Critical to the success of this implementation process is the commitment of the Laboratory Manager, Hematology Manager and the Trained Operator to provide the dedicated time, cooperation and coordination required to ensure completion of the necessary tasks.

## **Trained Operator Responsibilities**

- Provide setup information (such as reference intervals, decision rule criteria, critical limits, workflow details and LIS settings)
- Coordinate gathering of samples for reference intervals, method comparisons and truth tables
- Communicate with IT contact to ensure interface is ready as soon as possible
- Mentor your own laboratory staff
- Communicate with management to quickly address concerns and questions during process
- Complete required characterization studies
- Submit data for data analysis

The following is a synopsis of performance characteristics that are most studied. For more details, refer to Clinical Laboratory Standards Institute (CLSI).<sup>2</sup>

Some of the following procedures are completed one time at installation. Others will be required throughout the life of your instrument.

# Performed for Implementation and periodically (as per lab protocol and/or local regulatory agency)

- Verify Calibration using your instrument's commercially available calibrator
- Comparison studies against the previous analyzer or methods
- Manual differential comparisons
- Truth Tables for establishing and/or modifying laboratory flagging criteria
- Establishing Quality Control (QC) lab limits as per lab protocol or regulatory
- Verification of Adult Reference Intervals
- Verification of Body Fluid Analysis (if applicable)
- Verification of Retic results (if applicable)
- Verify accuracy and precision (Routine Quality Control procedures, including calibration verification)
- Verify Analytical Measuring Ranges (AMR)
- Verify comparison between automatic and manual modes (if applicable)
- Verify comparability between primary and back-up instrument

### Assistance from your Applications Specialist

The **characterization process** helps to assure that your new instrument, when used in your laboratory **by your testing personnel** for your patient population, **is performing as the manufacturer intended.** Your Applications Specialist will **assist** and **guide** you through the characterization steps; he or she will **not** perform the studies for you.

While your Applications Specialist is onsite, he or she (<u>with the assistance of the</u> <u>Trained Operator</u>) will complete the following:

- Calibration
- Setup of workstation for flagging limits and Decision Rules (if applicable) based on your laboratory protocols for review.
- Comparison study for CBC, automated differential and automated reticulocyte (if applicable). This study will compare your new instrument(s) to a single reference instrument.
- Verify the Analytical Measuring Range (AMR) using commercially available material per your laboratory protocols or regulatory agency.

All procedures will not be completed while your Applications Specialist is onsite. They may **guide you in the following**:

- Establishing Quality Control (QC) lab limits
- Verification of Body Fluid Analysis (if applicable)
- Verification of Retic results (if applicable)
- Truth Table analysis for comparison to manual differentials (if applicable)
- Verification of Adult Reference Interval
- Verify comparison between automatic and manual modes (if applicable)
- Verify comparability between primary and back-up instrument

In this Performance Verification Manual you will find discussions, procedures and worksheets to assist you with each stage of the characterization process. A Performance Verification Checklist is included to help you keep the process on track. Please note that the procedures provided outline the general steps for verifying the performance of your new system. Your current laboratory policies or your local regulatory agency may dictate more specific or additional rules for your laboratory to follow. We hope you find the information provided to be useful as you make the transition to your new hematology system.

Reagents received for the new system need to be removed from the shipping carton immediately and stored according to the storage criteria stated on the package insert located in the Beckman Coulter website. Some reagents are refrigerated and some are stored at room temperature.

<sup>1</sup> The Health Care Financing Administration (HCFA)'s State Operations Manual (Appendix C: survey Procedures and Interpretive Guidelines for Laboratories and Laboratory Services), published January 2016, details how to comply with the CLIA regulations.

Internet: http://www.hcfa.gov/pubforms/pub07/pub%5F07.htm

940 West Valley Road, Suite 1400 Wayne, PA 19087-1898 USA Phone: (610) 688-0100 Fax: (610) 688-0700 Internet: http://www.clsi.org/

<sup>&</sup>lt;sup>2</sup> CLSI

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFICATION MANUAL (Mm//dd/yy) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

#### Performance Verification CHECKLIST HEMATOLOGY COMPLETED DURING IMPLEMENTATION BY LABORATORY STAFF

Instru	nent/SN	Date Completed or N/A*	Tech	
	Hardware installation data verified			
	Precision verified			
	Accuracy verified			
	Carryover verified			
	Mode to Mode verified (if applicable)			
	Calibration			
	Set up DMS/Workstation/System Manager			
	Set up QC files			
	Enroll/Set up IQAP/eIQAP participant # for new instrument			
	Set up Interfacing to LIS (IT/LIS contact)			
	Setup Patient Flagging Limits and review criteria			
	Setup Decision Rules (if applicable)			
	Comparisons performed by laboratory staff			
	□ Samples run on new instrument and comparison method			
	□ CBC/Diff			
	Manual differentials for Truth Tables			
	Reticulocytes			
	Body Fluids			
	Comparison Data collated and submitted for Data Analysis			
	Verify Measuring Range (linearity)			
	Verify Reference Interval (Normal Ranges)			
	Establish QC lab limits (per lab protocol)			
	Verify Primary instrument to back up instrument (comparability)			
	Verify Specimen mixing study			
	Test Interfacing - Trained Operator and Lab IT Contact			
	Data analysis reports reviewed with appropriate Lab Staff			
	Pathology/Lab Director Sign off			
*0 .				

\*Some items non-applicable depending on Instrument, Test Menu, laboratory protocol and/or local regulatory agency

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFECATION MANUAL(mm/dd/yy) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

ASSESSMENT

## 2

# **Pre-Calibration Assessment**

#### Purpose

When performed under the following conditions:

• To assure the instrument is clean.

**PRE-CALIBRATION, CALIBRATION & QUALITY** 

- To assure the instrument is functioning properly prior to calibration.
- To assure when to verify vs. when to calibrate your instrument.
- To assure the instrument has the required reagents and enough of such to perform the procedure.

Summary and Explanation	This assessment should be done prior to running calibration. If this is the first time calibrating a new instrument, you should review how the commercial control is performing. If you are calibrating to verify calibration evaluate your current QC and the last 1-2 months of IQAP information for the calibrated parameters. Are you starting to trend high or low for any parameter? If you are high/low, how do you also compare to your peers? If your entire peer group is trending along with you this indicates that there may have been a change in the QC product itself. If you are trending, but not your peers, then it is time to address a possible calibration drift.
Principles of the Procedure	Pre-calibration assessment is the process of assuring that the instrument is in good working order and that the calibration factors before calibration are stable.

ProductReview and follow appropriate instrument Instruction for UseInformation(IFU) and calibrator package insert.

**Purpose** When performed under the following conditions:

- To assure that an instrument's data output accurately reflects sample input.
- To demonstrate that the performance of the instrument is consistent with the manufacturer's claims.
- To assure when to verify vs. when to calibrate your instrument.
- To bring an established method, device or analytical system online.
- To demonstrate acceptable performance as a follow-up to corrective actions taken after a failed proficiency-testing event.

#### Summary and Explanation

It is important that we remember that a Hematology calibrator cannot, by nature of the product, be an "absolute" constant. The product does contain a standard or exact number of cells at its creation and they are preserved to ensure a certain amount of time in which the product can be used. However, as these are viable cells they can undergo changes due to handling that can alter the final values. The parameters most commonly affected are the RBC (cells may lyse if frozen or overheated), MCV, and Plt (may increase if there is stroma from RBC lysis). It is good practice to evaluate your current QC.

	WBC	RBC	Hgb	HCT/MCV*	Plt	MPV
Level 1						
Level 2						
Level 3						
Cal Ref						
Cal mean						

\* Refer to specific commercial calibrator insert or instrument documentation for parameter to be calibrated.

Review each parameter

- indicate if it is running H/L
- decide which way you expect the calibration factor to change (increase or decrease)

**NOTE:** The higher the calibration factor the higher the final result. If you are trending on the high side, you would want the new cal factor to be smaller than the old factor.

If the new calibration factor does NOT make a change in the direction that you are expecting, even though it is for a parameter that you have determined needs to be recalibrated, do not change the factor at this time. Question the acceptability of the product. (Call Customer Support to troubleshoot)

Calibration Verification is done using commercial control. This allows you to demonstrate that the expected/needed changes are actually occurring to the QC material. If your instrument allows for entering comments to the QC file, document the Calibration by adding a comment to the last QC run before the calibration.

**Principles of** the Procedure The calibration procedure consists of comparing instrument measurements to known values for WBC, RBC, HGB, MCV, PLT and MPV. Calibration assures that an instrument's data output accurately reflects sample input. Calibration is performed using materials based on or traceable to known reference preparations or materials. In general, the procedure may indicate that the instrument requires standardization, by first determining the deviation from calibrator reference, and then applying recommended correction factors (CAL factors).

> The laboratory is responsible for the final calibration of the CBC parameters. Beckman Coulter recommends COULTER S-CAL calibrator, or an exact equivalent, as an acceptable alternative to whole blood calibration.

> In the normal process of tracking data for an extended period of time, your laboratory can make a specific decision to recalibrate a given parameter. Never adjust to a specific value for an individual sample.

> For best performance, verify and calibrate all the CBC parameters. The WBC differential, NRBC and Retic parameters are calibrated by an authorized Beckman Coulter Representative in your laboratory. The VCSn parameters do not require calibration in the laboratory.

**NOTE**: Ensure your SPM is properly maintained and the apertures are clean prior to calibration.

Refer to the calibration procedure listed in the instrument Instructions for Use (IFU) or on-line help.

#### When to Verify Calibration:

You should verify the calibration of your instrument:

- As dictated by your laboratory procedures, local or national regulations
- When controls begin to show evidence of unusual trends
- When controls exceed the manufacturer's defined acceptable limits
- If the average ambient room temperature changes more than 10°F or 12°C from the calibrating temperature.

If the procedure indicates you need to calibrate, continue with the calibration procedure.

#### When to Calibrate:

You should calibrate your instrument:

- At installation
- After the replacement of any component that involves dilution characteristics (such as the BSV) or the primary measurements (such as the apertures)
- When advised to do so by your Beckman Coulter Representative.
- If you fail verify calibration procedure.

ProductReview appropriate instrument Instruction for Use (IFU) and calibrator packageInformationinsert.

Specimen Collection and Preparation	Review appropriate instrument Instruction for Use (IFU) for requirements if applicable.
Materials Required But Not Provided	<ul> <li>Sufficient reagents to complete calibration.</li> <li>Product manuals, or Online HELP System.</li> <li>Appropriate Cell Control Kit.*</li> <li>Absorbent lint-free material such as paper wipes.</li> <li>Fresh, normal, whole-blood samples for pre-calibration procedures.*</li> <li>Calculator*</li> <li>Removable Media*</li> </ul> *Not required for all systems.
Procedure	Review appropriate instrument Instruction for Use (IFU) and calibrator package insert.
Calibration Details	Required calibration should be performed as stated in the calibrator package insert, and the instrument Instructions for Use (IFU).
Quality Control	After completing the calibration and verification procedures, the instrument is within the accuracy limit stated in your Product Manuals. Good laboratory practices recommend that each series of patient samples be preceded and followed by a quality-control check. Use COULTER cell controls to check the performance of the instrument. Verify the calibration with QC controls.
Results	The assigned values were established using representative samples from this lot of calibrator and are specific to the assay methodologies of the applicable Beckman Coulter hematology system reagents indicated in the calibrator Table of Expected Results (Assay Sheet).
	Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.
Limitations of the Procedure	Values assigned by other methodologies may be different. Such differences, if present, may be caused by inter-method bias.

Object Name: Hematology Performance Verification Manual (HPHEMATOLOGY PERFORMANCE EVERIFIC ATION<sup>8</sup>MANUAL Doc Title: Hematology Performance Verification Manual (HPVM)

Carryover	
Purpose	<ul> <li>When performed under the following conditions:</li> <li>To verify instrument performance.</li> <li>To verify instrument accuracy.</li> </ul>
Summary and Explanation	Carryover is verified by the installer on initial installation. A Carryover procedure may be a requirement of regulatory agency or part of laboratory protocol. This procedure may be used at any time to verify the performance of CBC, Diff and Retic parameters. <b>Most</b> Beckman Coulter hematology instruments have a <b>high to low</b> <b>carryover test</b> procedure that automates running samples and does the statistical calculations.
Principles of the Procedure	Carryover is the discrete amount of analyte carried by the measuring system from one specimen reaction into subsequent specimen reactions, thereby erroneously affecting the apparent amounts in subsequent specimen. Carryover is expressed conventionally as a percentage of the concentration of the analyte in the first specimen that is carried into the subsequent specimen.
Specimen Collection and Preparation	Review appropriate instrument Instruction for Use (IFU) for requirements if applicable.
Procedure	Review appropriate instrument Instruction for Use (IFU).
Results	When performing a carryover procedure, the calculated % carryover and/or background for each parameter is compared to the carryover and background limits for acceptability.
Limitations of the Procedure	Review appropriate instrument Instruction for Use (IFU) for requirements.

# Reproducibility/Repeatability

Purpose	When performed under the following conditions:
	<ul><li>To verify instrument performance.</li><li>To verify instrument precision.</li></ul>
Summary and Explanation	Reproducibility/Repeatability is verified by the installer on initial installation. A Reproducibility/Repeatability procedure may be a requirement of regulatory agency or part of laboratory protocol. This procedure may be used at any time to verify the precision of the instrument.
Principles of the Procedure	Reproducibility/Repeatability is a measure of the ability of the instrument to reproduce similar results when a sample is run repeatedly.
Specimen Collection and Preparation	Review appropriate instrument Instruction for Use (IFU) for requirements if applicable.
Procedure	Review appropriate instrument Instruction for Use (IFU).
Results	Reproducibility/Repeatability is assessed by replicate analysis of the same specimen $(n=10)$ . The closeness of agreement between the results of successive measurements of the same substance carried out under the same conditions of measurement.
Limitations of the Procedure	Review appropriate instrument Instruction for Use (IFU) for requirements.

## CBC PARAMETER COMPARISON 3

## Rationale

Method Comparison is a regulatory requirement to show agreement between two methods. Most often, a new instrument will be compared to the lab's current method. The comparison must be completed before the new method is placed into use. Patient specimens for analysis are usually acquired from the lab's routine population. Some specimens should be specially chosen to test the complete measuring range.

Before beginning the process of method comparison, verify the performance and calibration of both the current (reference) instrument and the test instrument. Both instruments should meet the manufacturer's specifications with regard to maintenance, background checks, reproducibility/repeatibility, carryover, calibration, and quality control. Refer to the product documentation for specifications.

## **Data Collection Guidelines**

The specimens collected for method comparison should reflect the typical patient population of the laboratory and be indicative of the clinically meaningful range for making medical decisions. This includes normal and abnormal inpatients, outpatients, pediatric and specialty clinics existing in the laboratory population. The quality of specimens collected is an important factor when collecting data. Poorly collected specimens contribute to poor data and lead to false conclusions. Properly collected and maintained specimens are:

- Collected in tubes containing a salt of EDTA and filled to the appropriate level to ensure the proper proportion of blood to anticoagulant.
- Collected according to the tube manufacturer's instructions.
- Have sufficient quantity to be run on both the current instrument and the test instrument.
- Processed on both instruments within one or two hours of each other and within the same time frame as the laboratory routinely analyzes specimens.
- Thoroughly mixed before processing to establish cellular equilibration.

# **General Considerations**

The number of specimens required for method comparison depends on the needs of the laboratory. CLSI (EP09)<sup>3</sup> recommends that at least 40 samples be analyzed. Samples chosen should span the upper and lower end of the measuring range using the available patient population.

The actual number of samples analyzed varies by laboratory and should be determined by the laboratory, considering the average number of samples run per day as well as the size of the laboratory. (i.e. a lab running 10,000 samples per day would not need to include 10,000 samples in the study.)

- Printed results for all specimens analyzed on all instruments should be obtained and labeled appropriately. For each specimen, collect the following data:
- Printout from the Evaluation instrument
- Printout from the reference instrumen
- Any confirmatory results
- If instrument has archive/export feature, archive/export specimen results to appropriate external media, such as a CD, 3.5 inch diskette or a flash drive and include with the printed results

# **Data Analysis**

Upon completion of specimen collection and processing, compare results from the current and test instruments and review for acceptable comparison. Methods of data analysis include:

- Calculation of mean difference statistics.
- Determination of method comparison and linear regression.
- Graphing.

J. Westgard and M. Hunt state "Statistical tests do not provide the criteria for acceptability. <u>Acceptability depends on whether the errors limit the clinical usefulness of the method.</u> Statistical tests can provide specific estimates of errors upon which judgements can be made, but they are not a substitute for judgments."<sup>2</sup>

#### IMPORTANT

The responsibility to determine acceptable performance of data lies with the evaluating laboratory.

Each laboratory is ultimately responsible for interpretation of comparison data and developing its own evaluation criteria for acceptability of results.

# DIFF PARAMETERS COMPARISON 4

## Rationale

Method Comparison is a regulatory requirement to show agreement between two methods. Most often a new instrument will be compared to the lab's current method. The comparison must be completed before the new method is placed into use. Patient specimens for analysis are usually acquired from the lab's routine population. Some specimens should be specifically chosen to test the complete measuring range.

Before beginning the process of method comparison, validate the performance of both the current (reference) instrument or manual differential and the test instrument. Both methodologies should meet the appropriate specifications for quality control. Refer to your appropriate laboratory protocols, or the respective product documentation.

## Automated vs Manual Diff Data Collection

When comparing the automated differential to the manual differential, ensure that the inherent variations of slide preparation are minimized by:

- Making quality smears.
- Staining with quality stain.
- Using optically clean microscopes.
- Having qualified technologists review the smear(s).
- Specimens are identified by lab number on both printouts and microscopic slides for future reference.

CLSI (H20)<sup>3</sup> recommends that two technologists each perform a 200 cell manual differential on two different slides (total of 400 cells analyzed by each technologist). More cells counted results in a more precise and accurate reference against which to judge the accuracy of the automated method. Automated differential systems analyze thousands of cells. The Rumke Binomial Distribution Table<sup>4</sup>, published by C.L. Rumke in 1978, illustrates that the statistical uncertainty of reference values is a direct function of the number of cells counted. The imprecision of the manual differential is especially pronounced with low numbers (e.g., monocyte, eosinophil and basophil percent).

- Printed results for all specimens analyzed on all instruments should be obtained and labeled appropriately. For each specimen, collect the following data:
- Printout from the Evaluation instrument
- Printout from the Reference instrument
- Any confirmatory results
- If instrument has archive/export feature, archive/export specimen results to appropriate external media, such as a CD, 3.5 inch diskette or a flash drive and include with the printed results

## **Data Analysis**

If a laboratory intends to compare a current automated differential to the new automated differential, then linear regression may be used. Keep the following in mind if you choose to do this:

- Compare similar technology (e.g. 5-part diff to 5-part diff or 3-part diff).
- Use the same samples on each instrument.
- Run the samples within 2 hours on both analyzers.
- "Many white cell count parameters including basophils and eosinophils often have correlation coefficients in the range of 0.2 to 0.5. Slopes and intercepts for such parameters are virtually worthless. The only important statistic resulting from this analysis is the bias at the upper and lower limits of the normal range." (EP Evaluator Report Interpretation Guide 8.0.0.165, David G. Rhoads Associates, Inc.)

Linear regression may not be an adequate statistical tool for analysis of the differential parameters due to the variability of the manual differential. Commonly used methods for determining clinically acceptable method comparison between automated differential and manual differentials are:

- Mean Difference
- Binomial Envelope Distribution

#### Mean Difference

The mean difference shows how much higher or lower the test instrument is compared to the current, or reference instrument results.

#### IMPORTANT

The responsibility to determine acceptable performance of data lies with the evaluating laboratory.
### **Binomial Envelope Distribution**

Binomial envelope distribution takes into account the imprecision of the manual differential. On the following pages are examples of binomial envelope graphs based on the Rumke Binomial Distribution Table<sup>4</sup>.

To plot the binomial envelope graph:

- 1. Select the plot that represents the number of cells counted on the manual differential vs. the automated differential.
  - a. Use the 0-100 scale for neutrophils (granulocytes) and lymphocytes.
  - b. Use the 0-30 scale for monocytes (mononuclears), eosinophils and basophils.
- 2. Each differential parameter should be plotted separately using the appropriate scale.
- 3. Plot the X-axis using either the manual diff value or the reference instrument value for a parameter.
- 4. Plot the Y-axis using the new method's value.

When the test method (new instrument differential) is compared to the reference method (manual differential or reference instrument), it may be considered accurate when approximately 95 % of the data is contained within the envelope area.

#### NOTE

On the monocyte (mononuclear) parameter plot, some points may fall outside the upper limit. Several publications have documented that monocyte counts derived from flow cytometry systems will be somewhat higher and more accurate than the routine 100-cell visual differential. Monocytes are particularly prone to poor distribution on the typical wedge smear.<sup>7</sup> Refer to the Commentary section of Manual versus Automated Differentials: Comments.

## **Rumke Binomial Distribution Table**

95 % Confidence Limits for various percentages of blood cells of a given type as determined by differential counts.

a	n = 100	n = 200	n = 500	n = 1000
0	0 - 2	0 - 2	0 - 1	0 - 1
1	0 - 6	0 - 4	0 - 3	0 - 2
2	0 - 8	0 - 6	0 - 4	1 - 4
3	0 - 9	1 - 7	1 - 5	2 - 5
4	1 - 10	1 - 8	2 - 7	2 - 6
5	1 - 12	2 - 10	3 - 8	3 - 7
6	2 - 13	3 - 11	4 - 9	4 - 8
7	2 - 14	3 - 12	4 - 10	5 - 9
8	3 - 16	4 - 13	5 - 11	6 - 10
9	4 - 17	5 - 14	6 - 12	7 - 11
10	4 - 18	6 - 16	7 - 13	8 - 13
15	8 - 24	10 - 21	11 - 19	12 - 18
20	12 - 30	14 - 27	16 - 24	17 - 23
25	16 - 35	19 - 32	21 - 30	22 - 28
30	21 - 40	23 - 37	26 - 35	27 - 33
35	25 - 46	28 - 43	30 - 40	32 - 39
40	30 - 51	33 - 48	35 - 45	36 - 44
45	35 - 56	37 - 53	40 - 50	41 - 49
50	39 - 61	42 - 58	45 - 55	46 - 54
60	49 - 70	52 - 67	55 - 65	56 - 64
70	60 - 79	63 - 77	65 - 74	67 - 73
75	65 - 84	68 - 81	70 - 79	72 - 78
80	70 - 88	73 - 86	76 - 84	77 - 83
90	82 - 96	84 - 94	87 - 93	87 - 92
100	96 - 100	98 - 100	99 - 100	99 - 100

**a** the observed percentage of cells of a given type **n** the total number of cells counted

# **DIFFERENTIAL BINOMIAL GRAPHS - Examples**

These are examples of typical graphs that are used to display and analyze your differential data. The graphs show the binomial envelopes used when comparing an automated differential to the manual differential or automated differential to automated differential. Use the appropriate envelope based on number of cells counted manually. When correlating an automated differential to an automated differential linear regression may be used.



Use this scale graph to plot Neutrophils (granulocytes) or Lymphocytes.



Use this scale graph to plot Monocytes or Mononuclears, Eosinophils, Basophils or NRBCs.

P/N: 4277076AC

## **DIFFERENTIAL LINEAR REGRESSION GRAPHS - Examples**

When method comparison using an automated differential to an automated differential linear regression may be used.



Use this scale graph to plot Neutrophils (granulocytes) or Lymphocytes.



Use this scale graph to plot Monocytes or Mononuclears, Eosinophils, Basophils or NRBCs.

# **RETIC PARAMETERS COMPARISON** 5

### Rationale

Method Comparison is a regulatory requirement to show agreement between two methods. Most often a new instrument will be compared to the lab's current method. The comparison must be completed before the new method is placed into use. Patient specimens for analysis are usually acquired from the lab's routine population. Some specimens should be specifically chosen to test the complete measuring range.

Before beginning the process of method comparison, verify the performance of both the current (reference) instrument or manual reticulocyte and the test instrument. Both methodologies should meet the appropriate specifications for quality control. Refer to your appropriate laboratory protocols or the respective product documentation for these procedures.

### **Reticulocyte Data Collection**

Collection of data should ensure that a typical distribution of normal patients and various abnormalities are included in the study. Obtain specimens with as wide a range of values as possible, over the reportable range of the method. Refer to your specific instrument documentation for ranges. Approximately one half of the specimens analyzed should be beyond the laboratory normal range if possible.<sup>28</sup>

The specimens should be:

- Collected according to the tube manufacturer's instructions.
- Analyzed on the automated system and prepared for reference counts within 2 hours of each other or as close in time to each other as possible. Analysis should occur within the same time frame as the laboratory routinely analyzes specimens.
- Stored according to instrument documentation for temperature and stability.

If using a reference manual reticulocyte count for comparison, ensure that the inherent variations of the method are minimized. These variations include but are not limited to:

- Preparation of stain
- Sample/stain mixture incubation
- Differences in slide preparation
- Area of review
- Technologist bias

CLSI (H44)<sup>28</sup> recommends two technologists each count 2000 RBCs using two different slides (total of 4000 RBCs counted). Data from each Beckman Coulter automated reticulocyte analysis are based on over 32,000 cells. More cells counted result in a more precise and accurate reference against which to judge the accuracy of the automated method. Statistical uncertainty of reference values is a direct function of the number of cells counted. The imprecision of the manual reticulocyte count is especially pronounced due to the low frequency of reticulocytes in the peripheral blood.

• Two counts should agree within the laboratory's acceptance limits. If the two counts do not agree, a third or "referee" count should be performed.

### **Reticulocyte Considerations**

- It is recommended that the determination is done promptly after collection of the blood specimen or alternatively, that the specimen be stored in such a way that it remains stable until the reticulocyte count is performed. With some reticulocyte dyes, there is an apparent *in vitro* maturation and subsequent disappearance of some of the reticulocytes, which is both time and temperature-dependent. For optimum performance, all specimens should be analyzed within time limits recommended by the manufacturer(s). The same specimen should be analyzed <u>within 2 hours</u> on all instruments for accurate comparison.
- If a manual retic count is performed, results of other test procedures, i.e. RBC morphology, should also be noted.

**NOTE:** Poor statistics could be obtained because of the imprecision of the manual reticulocyte count. This will be especially pronounced in cell populations with lower numbers.

- Specimens should be selected to test the clinical sensitivity of the instruments and to represent the typical distribution of normal and abnormal.
- 50 % of specimens representative of normal range specimens
- 50 % of specimens representative of abnormal specimens (25 % decreased and 25 % increased reticulocyte counts)
- Results should span as much of the clinical range of the instrument as possible.
- If flow cell clogs, incomplete computation, partial aspirations or results with System Messages are observed for any specimen, repeat the analysis of that specimen on the evaluation instrument. Submit both original and repeat print outs.
- Printed results for all specimens analyzed on all instruments should be obtained and labeled appropriately. For each specimen, collect the following data:
  - Printout from the Evaluation instrument
  - Printout from the Reference instrument
  - Manual reticulocyte results (if applicable)
  - Any confirmatory results
  - If instrument has archive/export feature, archive/export specimen results to appropriate external media, such as a CD, 3.5 inch diskette or a flash drive and include with the printed results

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFICATION MANUAL(*mm/dd/yy*) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

### **Data Analysis**

Upon completion of specimen collection and processing, compare results from the test instrument and current instrument or methodology and review for acceptable method comparison.

### IMPORTANT

Each laboratory is ultimately responsible for interpretation of method comparison data and develops its own evaluation criteria for acceptability of results.

Possible methods of data analysis include:

- Patient comparison, including:
  - Descriptive statistics (Mean, SD, CV)
  - Mean difference analysis (bias, accuracy)
  - ➤ Graphing

Linear regression, a commonly used statistical tool for laboratories, is not generally considered the best way to analyze Reticulocyte parameters where the manual reticulocyte count is used as reference. This is due to the high imprecision of the manual method and the low frequency of Reticulocytes in the peripheral blood. A more appropriate method of comparison is Mean difference analysis. Each individual laboratory must decide which method is appropriate for its own needs.

### TRUTH TABLE ANALYSIS 6

### Introduction

One of the most important tasks to accomplish during the Implementation of any automated differential analyzer is to establish/verify an effective flagging protocol. **The purpose of the flagging protocol is to identify those samples which require a slide review so that significant morphology detail can be added to the automated report**. The extent to which an analyzer can effectively screen "normal" vs "abnormal" is defined as the instrument's clinical sensitivity. The Truth Table is an effective tool for evaluation of the instrument's clinical sensitivity.

The first step in establishing an effective flagging protocol is to clearly define those findings that are considered clinically significant. A significant finding would be considered some detail that could add value to the report and ultimately affect/improve patient care.

The next step is to identify ways to flag for these findings. Some flags and codes are already built into the analyzer and are generated when the analyzer detects an unusual or unexpected distribution of particles. Additional flags can be set up by the operator to reflect the laboratory's review requirements and to support the screening process.

Specimens determined to be abnormal by the automated system should be evaluated according to the laboratory's protocols. Specimens determined to be normal by the automated system could be released without further review (autovalidation). The evaluation of the instrument's clinical sensitivity is critical because over-flagging may lead to unnecessary differential review, while underflagging may miss abnormal samples.

The manual differential is used as a reference method to assess whether the instrument's flagging is appropriate. The lab determines the flagging limits which separate normal results from abnormal results, classifying each sample as a True Negative, True Positive, False Negative or False Positive. Reviewing and possibly adjusting the operator definable flags may further improve the efficiency of the automated differential system.

#### **IMPORTANT NOTE**

Beckman Coulter, Inc. does not claim to identify every abnormality in all samples. Beckman Coulter, Inc. suggests using all available flagging options to optimize the sensitivity of instrument results based on your patient population. All flagging options include reference ranges (H/L), action and critical limits, definitive flags, suspect flags, system messages, parameter codes, delta checks, decision rules and system alarms. Beckman Coulter, Inc. recommends avoiding the use of single messages or outputs to summarize specimen results or patient conditions.

All Truth Tables and associated recommendations reflect the extent to which your workflow may be managed efficiently but are limited to the sample data submitted for evaluation. Accepting any of the suggested recommended changes to the current review criteria would require appropriate changes to action limits and/or lab protocol that were in place at the time of this evaluation. Finally, all observations are recommendations and subject to your review and discretion in the formation of your Laboratory's review and flagging protocols."

### Automated vs Manual Diff Data Collection

When comparing the automated differential to the manual differential, ensure that the inherent variations of slide preparation are minimized by:

- Making quality smears.
- Staining with quality stain.
- Using optically clean microscopes.
- Having qualified technologists review the smear(s).

Specimens are identified by lab number on both printouts and microscopic slides for future reference.

CLSI (H20)<sup>3</sup> recommends that two technologists each perform a 200 cell manual differential on two different slides (total of 400 cells analyzed). More cells counted results in a more precise and accurate reference against which to judge the accuracy of the automated method. Automated differential systems analyze thousands of cells. The following Rumke Binomial Distribution Table<sup>4</sup>, published by C.L. Rumke in 1978, illustrates that the statistical uncertainty of reference values is a direct function of the number of cells counted. The imprecision of the manual differential is especially pronounced with low numbers (e.g., monocyte, eosinophil and basophil percent).

# Truth Table Data Collection

- The ideal truth table would include a minimum of 100 samples with 50 % having an abnormal slide review. Specimens should be collected into K2 or K3 EDTA. For optimum performance, all specimens should be analyzed within time limits recommended by the manufacturer(s). Analyze the same specimen within 2 hour on all instruments for accurate comparison. Refer to CLSI Standard for Reference Leukocyte (WBC) Differential Count (Proportional) and Evaluation of Instrumental Methods; H20-A2 or current revision.
- In addition, a complete manual differential should be performed. A complete manual differential includes WBC differential, RBC and Plt morphology and WBC and Plt estimate. Results of other test procedures, e.g. manual Plt count, should be also noted.
- Poor statistics could be obtained because of the imprecision of the manual differential. This will be especially pronounced in cell populations with lower numbers.
- The statistical error of the manual differential can be reduced by performing several 200-cell differentials. Differential counts by several technologists are preferred. The CLSI protocol (H20) recommends two (2) technologists each doing a 200-cell differential for each sample. If you choose to have multiple technologists perform manual differentials, average them and submit the average differential counts.
- Ensure a random sampling of the population. Specimens analyzed should therefore represent your a general hospital population. Results should span as much of the clinical range of the instrument as possible. Check your instrument specifications.
- Random sampling should include:
  - > specimens with normal values.
  - abnormal specimens representing various types of leukocyte (WBC) disorders including but not limited to:
    - leukemias lymphocytosis lymphopenia
    - granulopenia granulocytosis
    - eosinophilia basophilia
  - abnormal specimens representing various types of erythrocyte (RBC) disorders including but not limited to:
    - hemoglobinopathy
    - polycythemia
    - microcytic anemias macrocytic anemias
  - abnormal specimens representing various types of platelet (Plt) disorders, including but not limited to:
    - thrombocytosis thrombocytopenia
    - morphological platelet disorders
- If vote-outs, incomplete computation, partial aspirations or results with System Messages are observed for any specimen, repeat the analysis of that specimen on the evaluation instrument. Submit both the original and repeated values.
- Printed results for all specimens analyzed on all instruments should be obtained and labeled appropriately. For each specimen, collect the following data:
  - Printout from the Evaluation instrument(s)
  - Printout from the reference instrument 0900324381e1d31a

- > Manual differential results, note the total number of cells counted
- > Any confirmatory results
- ➢ If instrument has archive feature, archive specimen results to appropriate external media, such as a floppy disk a CD, 3.5 inch diskette or a flash drive and include with the printed results

## **Define Review Limits and Review Protocol**

Purpose	To evaluate current or proposed review limits for your laboratory's patient population and specific needs.
Procedure	<ol> <li>Define and record a Test instrument review protocol on the Differential Review Limit Protocol Worksheet provided.</li> <li>Record the instrument-generated flags that require further action.</li> <li>Define the high and low limits your laboratory uses to take slide review action on an automated differential result.</li> <li>Be specific, define less than and greater than with "or equal to" if applicable; e.g. if &gt;10 % is your limit for Eosinophils, then 10 % is normal or negative and 11 % is abnormal or positive. Review appropriate instrument Instruction for Use (IFU).</li> </ol>
Results	Depending upon your protocol, results falling outside these abnormal limits require either a smear scan or a full manual differential.

# DIFFERENTIAL REVIEW LIMIT PROTOCOL WORKSHEET

- 1. Define Test Instrument Review Protocol (Abnormal Limits)
  - a. All instrument generated flags will be used. List exceptions below.
  - b. Operator Defined Flags

	Instrument Review Limits						
	Low	High					
WBC							
RBC							
Hgb							
НСТ							
MCV							
MCH							
MCHC							
RDW CV							
RDW SD							
Plt							
MPV							
Ne %							
Ly %							
Mo %							
Eo %							
Ba %							
NRBC %							
Ne #							
Ly #							
Mo #							
Eo #							
Ba #							
NRBC #							

2. Define Reference Manual Differential Abnormal Limits

Segmented Neutrophils	 Metamyelocytes	
Band Neutrophils	 Myelocytes	
Lymphocytes	 Promyelocytes	
Variant Lymphocytes	 Blasts	
Monocytes	 NRBCs	
Eosinophils	 RBC morphology	
Basophils	 Plt morphology	
	WBC morphology	

## How to establish a Truth Table

The automated differential's sensitivity, specificity and efficiency in distinguishing normal and abnormal specimens are determined using Truth Table analysis. The manual differential is used as a reference to assess whether the instrument has correctly classified the specimen as normal (released without further review-autovalidation) or abnormal (requiring review).

To perform a Truth Table Analysis:

- 1. For each specimen analyzed on the Beckman Coulter Analyzer:
  - a. Classify the instrument results as Normal ("negative") if no flags or messages are present.
  - b. Classify the instrument results as Abnormal ("positive") if flags or messages are present.
- 2. Perform a manual differential on all samples evaluated in Step 1.
  - a. Classify as Abnormal ("positive") any morphological or distributional abnormality observed.
  - b. Classify as Normal ("negative") any manual differential count in which all cell types are normal and within your established limits.
- 3. Categorize each specimen as one of the following four categories and record the results on the differential Truth Table Worksheet:
  - a. **True Negative (TN):** Normal (negative) by both test (new instrument) and reference (manual diff) methods.
  - b. **True Positive (TP):** Abnormal (positive) by both test (new instrument) and reference (manual diff) methods.
  - c. False Negative (FN): Normal (negative) by test method (new instrument) and abnormal (positive) by reference (manual diff) method.
  - d. **False Positive (FP):** Abnormal (positive) by test method (new instrument) and normal (negative) by reference method (manual diff).
- 4. Calculate and record on the Differential Truth Table Worksheet the following<sup>5</sup>:
  - a. **True Negative**: the percentage of specimens considered normal (negative) by both the test method (new instrument) and the reference method (manual diff).

### % TN = (# True Negatives / total number of specimens) x 100

b. **True Positive**: the percentage of specimens considered abnormal (positive) by both the test method (new instrument) and the reference method (manual diff).

### % TP = (# True Positives / total number of specimens) x 100

c. False Negative (FN): the percentage of specimens considered normal by the test method (new instrument) and abnormal by the reference method (manual diff).

### **Galen/Gambino**<sup>5</sup>

% FN = (# False Negatives / total number of specimens) x 100 CLSI

% FN= (True Negative/(True Positive + False Negative)) x 100

d. **False Positive**: the percentage of specimens considered abnormal by the test method (new instrument) and normal by the reference method (manual diff).

e.

Galen/Gambino<sup>5</sup>

% FP = (# False Positives / total number of specimens) x 100 CLSI

- % FP = (False Positive/(False Positive+True Negative)) X 100
- **Specificity**: the percentage of manual differential normals that were also normal on the instrument.

Specificity = # True Negatives ÷ # (True Negative+False Positives) x 100

f. Sensitivity: the percentage of manual differential abnormals that were also abnormal on the instrument.
 Sensitivity = # True Positives ÷ # (True Positives+False)

Negatives) x 100

g. **Predictive Value of a Negative Test (PVN):** the percentage of specimens that were normal on both the instrument and the manual differential.

PVN = # True Negatives ÷ # (True Negatives + False Negatives) x 100

h. **Predictive Value of a Positive Test (PVP)**: the percentage of specimens that were abnormal on both the instrument and the manual differential.

**PVP** = # True Positives ÷ # (True Positives + False Positives) x 100

i. **Agreement**: the percentage of specimens correctly categorized by the automated differential.

[# (True Positives + True Negatives) ÷ total number of specimens] x 100

### IMPORTANT

The Truth Table format can be used for a variety of purposes. It is important to understand that the statistics reflect the patient mix that is used during the study. The focus of an Implementation Truth Table is to optimize Sensitivity. Focus on the False Negative Samples. Evaluate if there is a particular flag that could be added or modified so that these same samples would be flagged. Determine what is an acceptable balance between False Negative (missed significant findings) and False Positive samples (increased review rate). The responsibility to determine acceptable performance of data lies with the evaluating laboratory.

# DIFFERENTIAL TRUTH TABLE WORKSHEET

1. Classify differentials into one of four categories and total each column.

	Reference (Manual Differential)							
Tost		Normal (Negative)	Abnormal (Positive)	Total				
(New Instrument)	Normal (Negative)	(TN)	(FN)					
	Abnormal (Positive)	(FP)	(TP)					
	Total							

2. Calculate the following parameters.

PARAMETER	CALCULATION	RESULT (%)
% TN	Number of TN/Total	
% TP	Number of TP/Total	
% FN	Number of FN/Total	
% FP	Number of FP/Total	
Specificity	[# TN/ # (TN + FP)] x 100	
Sensitivity	[# TP/ # (TP + FN)] x 100	
Predictive Value of a Negative Test	[# TN/ # (TN + FN)] x 100	
Predictive Value of a Positive Test	[# TP/ # (TP + FP)] x 100	
Agreement	[# TP + # TN/Total] x 100	

3. Reviewed by \_\_\_\_\_ Date \_\_\_\_\_

## MANUAL DIFF WORKSHEET

Sample ID		Average	Your lab's definition of Positive**
Tech Initials	Tech Initials		
Segs	Segs	Segs	Segs
Bands	Bands	Bands	Bands
Meta	Meta	Meta	Meta
Myelo	Myelo	Myelo	Myelo
Pro	Pro	Pro	Pro
Blast	Blast	Blast	Blast
Lymph	Lymph	Lymph	Lymph
Atyp Ly	Atyp Ly	Atyp Ly	Atyp Ly
Mono	Mono	Mono	Mono
Eos	Eos	Eos	Eos
Baso	Baso	Baso	Baso
NRBC	NRBC	NRBC	NRBC
Other	Other	Other	Other
Total	Total	Total	
Significant RBC morph	Significant RBC morph	Significant RBC morph	Significant RBC morph
Significant PLT observations	Significant PLT observations	observations	Significant PLT observations
		Determine if your sample is	"Positive" or "Negative"
		by comparing your final (av	erage) result with your lab's
		definition of "positive" (act	ually seen on slide)
**Input lab's limits for a master		Circle one	Positive Negative
and copy as needed.		$  \rightarrow$	

Object Name: Hematology Performance Verification Manual (HPVM)	Effective:	09/17/18 (mm/dd/yy)
Dochen Andrew Bertermance Wentean on Manual (HPWM) ANUAL	Status:	Released

## OTHER STUDIES 7

This section contains guidelines and examples. Each lab must establish their own protocols for instrument implementation and use. The following studies may be required by your current laboratory policies or your local regulatory agency.

- **Body Fluids Comparison**
- □ Measuring Range (Linearity)
- □ Mode to Mode (if instrument has separate aspiration pathways)
- **Establish QC Lab Limits**
- **Gamma Reference Interval (Normal Range)**
- **General Science Study**

The following procedures will assist you in the evaluation a new instrument or clinical test. Not all of these procedures are required by CLIA; however, some of them may be required by the College of American Pathologists (CAP), the Joint Commission of Accreditation of Healthcare Organizations (known as JCAHO), and/or individual agencies. It is important to know local requirements as well as those of any other accrediting agencies that will impact your laboratory.

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFICATION MANUAL (mm/dd/yy) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

# **BODY FLUID Comparison**

Purpose	Method Comparison is a regulatory requirement to show agreement between two methods. Most often a new instrument will be compared to the lab's current method. The comparison must be completed before the new method is placed into use.						
Summary and Explanation	Beginning the process of method comparison, verify the performance of both the current (reference) instrument or manual cell count and the test instrument. Both methodologies should meet the appropriate specifications for quality control. Refer to your appropriate laboratory protocols or the respective product documentation for these procedures.						
Principles of the Procedure	A corrected RBC (RBC minus WBC) is often necessary in order to achieve acceptable method comparison. The RBC count will include all particles counted that are greater than 36 fL. When the WBC is elevated to a point of statistical relevance correct the RBC for the WBC count.						
Specimen Collection and Preparation	Patient specimens for analysis are usually acquired from the lab's routine population. Some specimens should be specifically chosen to test the complete measuring range.						
Results	Upon completion of specimen collection and processing, compare results from the current and test instrument or methodology and review for acceptable method comparison.						
	<b>IMPORTANT</b> Each laboratory is ultimately responsible for interpretation of method comparison data and develops its own evaluation criteria for acceptability of results.						
	Possible methods of data analysis include:						
	<ul> <li>Patient comparison, including:</li> <li>Descriptive statistics (Mean, SD, CV)</li> <li>Mean difference analysis (bias, accuracy)</li> <li>Graphing</li> </ul>						
Limitations of the Procedure	Linear regression, a commonly used statistical tool for laboratories, is not generally considered the best way to analyze parameters where the manual count is used as reference. This is due to the high imprecision of the manual method and the low frequency of cells in some fluids.						

# **Body Fluid Data Collection**

			Manual		Automated count			Background	
#	Sample ID	Fluid Type	WBC	RBC	TNC	RBC	corrected RBC	TNC	RBC
1									
2									
3									
4									
5									
6									
7									
8									
9									
10									
11									
12									
13									
14									
15									
16									
17									
18									
19									
20									
21									
22									
23									
24									
25									
26									
27									
28									
29									
30									

Object Name: Hematology Performance Verification Manual (HPVM)	Effective:	09/17/18 (mm/dd/yy)
DEFME PENAGISEPERORMANCE VERICARIAN AND AL	Status:	Released

## Measuring Range (LINEARITY)

Purpose	Verification own laborato use.	of linearity ma ry protocol. Th	by be required by your regulatory agency or your his procedure verifies instrument specific intended						
Summary and Explanation	You may purchase a <b>commercial product</b> . Follow the product insert for use. If you do not purchase a commercial product you may use this alternate method using whole blood, making serial dilutions of selected samples. Whole blood linearity testing involves analysis of diluted specimens to measure performance of WBC, RBC, HGB, and PLT throughout the linear range. The specimens required will be instrument specific based on the manufacturer								
Materials Required But Not Provided	<ul> <li>Comi</li> <li>Whol</li> <li>Clear</li> <li>Isotor</li> <li>Gradu</li> </ul>	Commercial product Whole Blood Clean Tubes Isotonic diluent or blood bank saline Graduated pipette							
Procedure	The followin linearity: 1. For each p the linear preparation 2. A concent > 60-65 % 3. Label five 4. Use isotor should be 5. Prepare di 6. Ensure the the prepara Prepare the s <b>Dilution</b> 100 %	ng steps are a barameter, select range. Ensure th n and aspiration rated specimen 6. clean tubes 80 nic diluent or bl run as soon as lutions using th e 100 % sample ation process. amples as follo <b>Specimen</b> 10 parts	general procedure that may be used to verify et an appropriate specimen at/near the high end of hat sufficient specimen is collected for dilution n. a can be used. Specimens should not have a Het 0 %, 60 %, 40 %, 20 % and 0. lood bank saline to make dilutions. Samples possible after preparation. he original specimen as 100 %. e and the dilutions remain well mixed throughout ows: Diluent 0 parts						
	80 % 60 % 40 % 20 % 10 % 0 %	8 parts 6 parts 4 parts 2 parts 1 part 0 parts	2 parts 4 parts 6 parts 8 parts 9 parts 10 parts						
	7. Analyze each of the dilutions in triplicate.								

8. Record the results and plot on XY graph.

**Results** Visually examine the XY plot. The data should appear linear with no outlying points.

#### Tips:

- You may use your instrument's Reproducibility/Repeatability screen, CBC only mode
- If your instrument uses blood detectors, remember that the highest sample dilutions (0, 10, 20 %) may cause aspiration errors because the Hemoglobin is below the lower threshold on some instrument models. If you have the option, disable blood detectors for these dilutions so that the samples may be run in the automatic or primary test mode. Refer to specific instrument model documentation.
- Master Worksheets are located in Section 11, you may make copies.

### **Example Graph for WBC**



LimitationsCheck your instrument documentation for specifications to verify its linearity<br/>limits.Procedure

## MODE TO MODE COMPARISON

Purpose	This procedure is used on instruments that have different aspiration pathways. The procedure is used to verify that there are no significant differences between the manual and automatic aspiration modes, ensuring the WBC, RBC, Plt, and Hgb results for a specimen are the same regardless of the mode selected.					
Summary and Explanation	Beckman Coulter recommends that Quality Control checks be performed using patient or commercial controls in both automatic (primary) and manual (secondary) modes at intervals established by your lab.					
Product Information	Refer to your own laboratory protocol or your local regulatory agency. When using commercial control, refer to the package insert to determine which mode to use.					
Specimen Collection and Preparation	Collect the necessary data in a Reproducibility mode if available on your instrument; run samples in the manual mode and again in the automatic mode. If using a Reproducibility mode ensure you print the runs before deleting for the second run set.					
Materials Required But Not Provided	10 Fresh normal blood specimens					
Procedure	<ol> <li>Select your instruments CBC mode.</li> <li>Cycle the samples, aspirating each sample once.</li> <li>After analyzing all the samples, obtain the mean values in both aspiration modes. If running in Reproducibility the means will be calculated.</li> </ol>					
Results	Compare the results to the specifications in your product documentation.					
	Calculate the absolute and percent difference using these formulas:					
	Absolute Difference = Automatic Mode Mean minus Manual Mode Mean					
	Percent Difference = (Absolute Difference / Automatic Mode Mean) x 100					
	If the % and absolute difference (whichever is greater) results are within tolerance, you have verified the current mode to mode calibration is correct.					
Limitations of the Procedure	If results do not meet your instrument specifications contact Customer Support.					

## MODE TO MODE CONSIDERATIONS

The mode to mode specifications stated in your instrument documentation state a range for absolute difference and percent difference (use whichever is greater). By using the procedure stated and running only 10 samples each, you are either verifying that no differences exist, or identifying the need for further investigation. If required, the actual calibration procedures require 50 samples run in both modes or 10 samples run in triplicate.

It is important you have clarity from your regulatory agency. Ask for written guidelines, for example:

CAP Hematology - Coagulation Checklist: 07/28/2015 Sampling Mode Comparison HEM. 30070 Phase I YES NO

There are records that at least annually compare all results obtained for patient specimens analyzed in the multiple sampling modes of the CBC analyzer (e.g. "open" and "closed" modes) to ensure that they are in agreement?

#### NOTE

Different modes may involve a different sample path before analysis. When samples are analyzed in more than one mode, it is important to ensure that all modes function properly. Reanalysis of a previously analyzed sample should be performed in the alternate mode(s), and results should agree with the initial mode within the tolerance limits established for agreement by the hematology laboratory's quality control program, and any recommendations by the instrument manufacturer. Mode-to-mode correlation is not necessary for those analyzers which use the same pathway for all modes.

#### **Evidence of Compliance:**

 $\ensuremath{\mathbbmath$\mathbbms$}$  Written procedure for sampling mode comparison with defined criteria for agreement AND

2 Records of sampling mode comparison studies

### **Quality Control Test - Hematology Analyzers with Two Sample Modes**

http://www.jointcommission.org/AccreditationPrograms/LaboratoryServices/Standards/09 FAQ

Updated | April 11, 2016

**Q:** What is required for daily QC in hematology for an analyzer with two sample modes?

A: Since there are two distinct sample pathways, QC is required for each sample mode according to the parameters established in the hematology standards.

- Performance of at least two levels of commercial controls every 24-hours of patient testing
- The controls used in the 24-hour period test the entire range of reported results
- Patient controls may be used to supplement the commercial controls if an acceptable level of precision has been defined.

For a single or primary mode, must assay at least two levels of quality control material each twenty-four hours of patient testing. In addition, the combination of controls must span the reportable patient range. Many laboratories still choose to perform three levels of commercial controls, with at least one level run on each 8-hour shift. In lieu of using commercial controls for QC of the secondary mode, laboratories may instead use patient controls whose values have been established in the primary mode when commercial controls were within acceptable ranges. For satisfactory compliance, patient controls should test the reportable range and have acceptable levels of precision defined in policy.

### MODE TO MODE COMPARISON LOG

Instrument:	Serial Number:

Acceptable Limits\*\* WBC \_\_\_\_ RBC \_\_\_\_ HBG \_\_\_\_ MCV \_\_\_\_ PLT \_\_\_\_

\*\*Lab can define their own acceptable lab limits. Reference the mode to mode specifications in your instrument documentation for manufacturer claims.

Date:	Aspiration Mode	WBC	RBC	HGB	MCV	PLT
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference					
	Automatic					
	Manual					
	difference	1				
	Automatic					
	Manual	1	1			
	difference					
# ESTABLISH LABORATORY LIMITS FOR QC

## Purpose of Establishing Laboratory QC Limits

As noted by Westgard, <sup>24</sup> "Means, standard deviations, ranges and other data from outside your laboratory does not reflect the individual, particular conditions of your lab. The use of data supplied from outside the laboratory to provide means, standard deviations and control ranges is meant to be a temporary workaround." Thus each laboratory must establish its own commercial control means and ranges, using a cumulative approach to calculations.

It is important you have clarity from your regulatory agency. Ask for written guidelines. The following procedure is suggested by JCAHO.

- 1. Obtain your last ten IQAP reports or your cumulative control statistics for the past 10-12 months printed from your instrument.
- Document your CV % for each directly measured parameter for each level of control. (Your IQAP report lists CV % for directly measured parameters only. Other parameters are considered "calculated" parameters). See example below:

## CAUTION

CV % from IQAP reports or cumulative statistics printouts may be erroneously low if your lab routinely deletes out-of-range controls.

Month #	Abnormal II Hgb CV %	Abnormal I Hgb CV %	Normal Hgb CV %
1	1.29	0.80	0.41
2	1.17	0.69	0.52
3	1.23	0.71	0.49
4	1.25	0.66	0.70
5	0.98	0.59	0.65
6	1.32	0.48	0.71
7	1.27	0.62	0.51
8	1.38	0.74	0.79
9	1.30	0.68	0.48
10	1.15	0.79	0.62
Average	1.234	0.676	0.588

## EXAMPLE

- 3. Compute the average CV %
- 4. From the control package insert, locate the Assigned Value for each parameter.

5. Multiply the average CV % obtained in step 3 above by the Assigned Value then divide by 100. This is 1 SD. The formula is: (CV % x Assigned Value) / 100 = 1 SD.

## EXAMPLE

Month #	Abnormal II	Abnormal I	Normal	
	Hgb CV %	Hgb CV %	Hgb CV %	
1	1.29	0.80	0.41	
2	1.17	0.69	0.52	
3	1.23	0.71	0.49	
4	1.25	0.66	0.70	
5	0.98	0.59	0.65	
6	1.32	0.48	0.71	
7	1.27	0.62	0.51	
8	1.38	0.74	0.79	
9	1.30	0.68	0.48	
10	1.15	0.79	0.62	
Average CV %	1.234	0.676	0.588	
Assigned Value	6.7	12.3	17.0	
1 SD (CV % x Assigned Value) / 100	0.083	0.083	0.10	
3 SD (1 SD x 3)	0.25	0.25	0.30	

6. Multiply 1 SD by 3 to obtain 3 SD. Enter the 3 SD as the expected range (Lab Limits) in the control setup screen. Enter package insert expected ranges for "calculated" parameters.

### NOTE

• VERY IMPORTANT - Evaluate the 3 SD ranges derived from this procedure! If the 3 SD ranges are too narrow and is not clinically relevant, the range may be widened with the authorization of the laboratory (medical) director to a range more reflective of the clinical applications of the parameter. The derived 3 SD range should not exceed the manufacturer's package insert range or the HCFA allowable variation for the parameter. Retain documentation for any changes authorized by the director.

• The manufacturer's package insert ranges may be used if there is verification that the mean obtained by the laboratory reflects the manufacturer's mean AND if the laboratory director assures that the manufacturer's range is narrow enough to provide results with meaningful clinical applications.

• There is no need to re-compute the average CV % unless changes in instrument precision or the control product have been noted. 0900324381e1d31a

## **ESTABLISHING LAB LIMITS CONSIDERATIONS**

The Joint Commission Use of Manufacturer Quality Control Ranges

http://www.jointcommission.org

Updated November 08, 2017

# Q; Can the manufacturer's stated ranges be used as the laboratory's quality control ranges?

A: The standards require each laboratory to establish their own control ranges through repetitive testing. However, there is an allowance to use manufacturer ranges when the following conditions are met:

- 1. The stated values correspond to the method and instrument used by the laboratory, and
- 2. The mean obtained by the laboratory reflects the manufacturer's stated mean, and
- 3. The Laboratory Medical Director assures the range is narrow enough to detect clinically significant error.

Manufacturer ranges may also be implemented if a test is used so infrequently that calculation of valid statistics is not possible. In settings where there is a high reproducibility (precise instrumentation, limited testing personnel), the laboratory's own calculated standard deviation (SD) may be small. When compared with the manufacturer ranges, a laboratory may find that the range spans more than the commonly used + / - 2 SD. Using the laboratory's calculated + / - 2 SD may produce unnecessarily narrow ranges, causing the testing personnel to frequently repeat QC and investigate when the controls performs outside laboratory's range, but within the manufacturer's range. Alternately, the full manufacturer range may be too broad to promote the detection of clinically significant error. Selection of the appropriate range is a balance between these two ends of the spectrum.

It is at the determination of the Laboratory Medical Director to approve quality control ranges after giving consideration to the clinically significant variance as compared to the statistically derived SD.

## CAP Requirements: Hematology - Coagulation Checklist: 07/28/2015

## HEM.20035 Phase II

Are tolerance limits (numeric and/or non-numeric) fully defined and documented for all hematology and coagulation control procedures?

**NOTE:** The goal is to have scientifically valid, logical "action limits" for quality control procedures that promptly alert the technologist of the need for immediate evaluation of the particular assay, including initiation of corrective action, before release of patient results.

## HEM.25870 Phase II

If commercially ASSAYED controls are used for CBC instruments, do control values correspond to the methodology and have target values (mean and QC ranges) are verified or established by the laboratory?

**NOTE:** Most commercial controls have expected recovery ranges for each parameter, provided by the manufacturer. The mean of such ranges may not be the exact target value in a given laboratory. Each laboratory should assign its own initial target value, based on initial analysis of the material; this target value should fall within the recovery range supplied by the manufacturer, but need not exactly match the package insert mean. The laboratory should establish specific recovery ranges that accommodate known changes in product attributes, assuming that calibration status has not changed.

**Other References:** CLSI Document C24 Internal Quality controls Testing: Principles and Definitions for establishing each laboratory's individual mean target values and ranges<sup>32</sup>, and CLIA Title 42 part 493 subpart K, Quality Control Test for Moderate and High Complexity, Section 493.1256, Standard Control Procedures.

# **REFERENCE INTERVALS (NORMAL RANGE)**

The CLIA Quality System Regulations became effective on April 24, 2003. Now the laboratory is required to check (verify) the manufacturer's performance specifications provided in the package insert – for accuracy, precision, reportable range, and **reference range**.

In order to **establish** a reference interval/range, Clinical and Laboratory Standards Institute (CLSI) recommends labs collect a sufficient number of **qualified** reference donors to produce a minimum of 120 samples. This can be a difficult task. Therefore, CLSI strongly encourages laboratories to **verify** reference values established by the manufacturer or established by another laboratory.<sup>28</sup>

To **verify** a reference interval/range:

- 1. If a lab has previously established a reference interval (using 120 samples) for its population, it may verify that reference interval by **transference**.
  - a. You may use transference if the original study was conducted using the same analytical system/reagents and the original study was conducted with similar subject population as your lab.
  - b. Perform a method comparison between your new method and the method used in the original study. You may use samples that exceed the reference interval. Evaluate the relationship between the two methods. If the new method provides values that are acceptably comparable, then the reference interval/range can be transferred.
- 2. Or, a lab may verify a reference interval/range established elsewhere or by the manufacturer by collecting a minimum of 20 samples from qualified, healthy reference donors. Your 20 donors must fairly represent your population and the population in the original study.

The original study's reference interval/range may be considered valid if no more than two of the 20 donors' values (or 10 % of the results) fall outside of the original study's limits.

In general, specimens should be handled in a manner similar to the patient samples.

**TIP**: You may use Control Folders to assist in gathering and computing reference range data. Verify the number of runs that a control folder will hold on your instrument. You may verify all adult normals or each sex individually. Set up a control folder for "Male" and another control file for "Female" or one folder called "Adult" for both sexes. It is not necessary to enter Assigned Values or Expected Ranges. Your instrument may have a "Studies" mode to use instead of Control Folders.

Carefully choose your samples to be included in the Reference Interval Study. If you have set up control folders on your instrument run them into the appropriate control folders. If you are using "Studies", you will also need to select batch. Ensure the samples are batched properly or run into correct control folders for male, female or adult.

Once you have **20** samples in the control folder(s), archive or export the control folders into a spreadsheet. Apply 95 % confidence limits to the data.

Example of formulas if you are using an Excel spreadsheet:

Lower 95 % confidence limits: = PERCENTILE (A2:A41, 0.025) Upper 95 % confidence limits: = PERCENTILE (A2:A41, 0.975)

## NOTE

A2 and A41 are examples of a range of data cells in the spreadsheet. Your actual range of data cells may be different.

It is up to your laboratory to evaluate the Reference Range data and determine the acceptability and utility of the data.

#### References

1. Clinical and Laboratory Standards Institute. Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline. CLSI document EP28-A3C, 2010

2. Clinical and Laboratory Standards Institute. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline. EP09-A3, 2013. 0900324381e1d31a

# **SPECIMEN MIXING**

**Purpose** The accuracy of whole blood cell counting by any method is based on the assumption that the specimen being analyzed is homogeneous. It is therefore, imperative that one ensures that a specimen is properly mixed prior to analysis. Each laboratory should have defined protocols for all pre-analytical (specimen collection, storage and mixing) and analytical (instrument) processes. These protocols should take into account the manufacturer's instructions for use for both the blood collection device and the analytical instrument.

Summary and Explanation When an anti-coagulated whole blood specimen tube is placed in a vertical position, the cellular components tend to fall to the bottom of the container. The erythrocytes (RBCs) concentrate at the bottom of the tube with leukocytes (WBCs) and platelets (PLTs) layered on the top of the cellular component. The WBC and PLT layer is commonly called the buffy coat. Plasma and cellular factors can affect the rate of settling. Elevated levels of fibrinogen, globulins and cholesterol have been reported to accelerate the rate of sedimentation. Anemia (low number of RBCs) also increases the rate of sedimentation of a sample. Specimens with high concentration of RBCs will settle slowly. The longer that the specimen has been allowed to settle, the longer the required time to resuspend the cells thoroughly.

> Adequate sample mixing is a critical part of an accurate cell counting method. A laboratory is best served to define their typical pre-analytical workflow and to verify that procedures are in place to ensure that samples received for cell counting are adequately mixed prior to analysis.

## Principles of the Procedure

If a settled specimen is not adequately mixed, results will be compromised. The pattern of results is dependent on how long the tube has settled and where the sample is taken from inside the tube. If the aspiration device is inserted through the buffy coat towards the bottom of the tube, there is a typical pattern of parameter results, which show elevated RBC/HGB and decreased WBC/PLT. See the following examples.

	WBC	RBC	HGB	HCT	MCV	MCH	MCHC	RDW	PLT	MPV
Example 1										
Run 1 Unmixed Specimen	7.0	5.88	16.6	50.9	86.5	28.3	32.7	20.5	171	7.4
Run 2 Actual results <sup>A</sup>	10.1	2.33	6.8	20.4	87.8 <sup>B</sup>	29.3	33.4	21.0	453	8.3 <sup>C</sup>
Relationship of Run 1 to Run 2	Û	Û	Û	Û	⇔	⇔	\$	\$	Û	⇔



A: Specimen in both runs remixed, 16 hours later

 $\square$  Decreased,  $\square$  Increased, No change [mixed and unmixed results essentially equal]

**B:** Slight increase due to age of specimen

C: Slight increase due to age of specimen

	WBC	RBC	HGB	MCV	PLT
Example 2					
Poorly mixed	1.5	2.43	8.0	92.3	8.0
Repeat after Mixing	1.8	1.48	4.7	92.1	10.0
Example 3					
Poorly mixed	5.6	5.64	17.0	86.1	192
Repeat after mixing	7.0	4.05	12.2	86.5	321
Example 4					
Poorly mixed	1.2	7.67	23.4	92.0	34
Repeat after mixing	3.5	4.15	12.8	92.5	171

If the sample is taken from the top of a settled tube, you may see the inverse effect with falsely high WBC/PLT and falsely low RBC/HGB. The effect on the percentage parameters such as the WBC Differential is unpredictable but typically not affected.

A mixing error may not be recognized immediately because the result set may be credible. This type of error is recognized if the laboratory performs delta checks by comparing the current result to a previous result. The delta review may be done by the instrument data management system, the laboratory information system or manual review of patient history. When the error has been acknowledged, one might suspect that there has been a sample mix-up because the results are credible. By reviewing the results of the parameters that are not count related, you will typically see that the MCV, MCH, MCHC and MPV correlate well between the runs. As noted in Example 1 above, the sizing parameters [MCV and MPV] may show slight increases as the sample ages. Sample results on the original tube may be compromised depending on the volume of blood that has been aspirated from the sample tube and the degree of mixing that was performed. The best strategy for verifying the results would be to redraw and rerun the sample.

Factors other than the blood specimen also contribute to inadequate mixing. Current venous blood drawing tubes utilize internal vacuum to draw the volume of blood proportional to the anti-coagulant concentration. The vacuum allows for an air space between the blood and tube stopper. Upon inversion, the air bubble helps the blood mix. If the tube stopper is removed and the tube is manually filled, the air space may not be adequate to mix the sample properly no matter how long it is mixed.

The most effective means to properly mix a whole blood specimen is by repeated inversion of the tube. A rotary type mixer, shaped like a wheel, can accomplish this mixing without intervention. Some tilt-type mixers can cause problems if the rocking speed is too slow or if a tube is placed at the tilt axis instead of the edge of the rocker bed. The best way to ensure that a specimeneis<sub>31a</sub>

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFECTION MANUAL (mm/dd/yy) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

adequately mixed is to invert the sample tube manually and observe the bottom and sides of the tube. If cells remain adhered to the inside of the tube, more mixing is necessary.

Automated Mixing Systems Automatic specimen tube processing for Coulter hematology systems was introduced in 1987 with the COULTER® STKR<sup>TM</sup> Hematology Analyzer. Specimen tubes in a cassette are placed in the loading bay of the instrument. These cassettes are then automatically transported, mixed, aspirated and analyzed. This same type of sample transport has been implemented on other Coulter analyzers such as the COULTER STKS, ONYX Autoloader, MAXM Autoloader, GEN•S, HMX Autoloader, LH Series, and DxH Series.

> As part of each new Beckman Coulter product validation, extensive studies are done internally to verify that the mixing of the specimen on the Coulter Hematology analyzer is sufficient to ensure that samples are homogeneous at the time of sample aspiration.

Studies may include the following types of challenges:

- Centrifuging specimens to simulate exaggerated settling times
- Concentrating specimens to simulate high RBC counts by removing plasma from the tube
- Various tube sizes and shapes (surface area to volume ratios)
- Different cap shapes and materials (standard latex, Hemogard type)
- Different anti-coagulants [K2, K3 and NA2EDTA] and concentrations
- Varying the volume of blood within a sample tube (over and under anticoagulation)
- Different vendor sample tubes
- Plastic versus glass
- Different storage temperatures

The efficacy of mixing is tested against the mode to mode specification for WBC, RBC, HGB and PLT and the accuracy specification for MCV, RDW and MPV. The test case results are compared to specimens mixed 5-10 minutes on a Clay Adams Nutator mixer prior to analysis on the hematology analyzer.

All cassettes in these studies were introduced and completed the entire premixing prior to processing. Even if properly processed, some tubes exhibit mixing problems. If additional pre-mixing is necessary for a particular tube type, it is stated in the Tube List in the Operators Guide. If samples are introduced directly, sample mixing may be inadequate to ensure quality results. Method No matter how comprehensive, internal testing cannot mimic all possible laboratory workflows. For that reason, a laboratory should perform a study with its own typical pre-analysis conditions. When setting up a study, considerations should include but may not be limited to:

- Specimen collection tubes types
- Tube size and material (glass versus plastic)
- Mode of collection (venous, syringe fill, micro-collection)
- Volume of blood collected within the sample tube
- Time between draw and analysis
- Storage temperature
- Transport mode to laboratory
- Mode of introduction of cassette to analyzer

Laboratory accreditation agencies look for quality assurance documentation for blood collection from point of patient identification to reporting results. Specimen mixing is just one of the pre-analytical steps in the process that needs to be documented as part of a complete quality assurance program, the CAP Accreditation Program Checklist<sup>4</sup>.

## **Commentary** Question: Heme:22000 [Phase II] asks:

Are all blood specimens collected in anticoagulant for hematology testing mixed thoroughly immediately before analysis?

They also include the following commentary.

Specimens collected in anticoagulants for hematology studies must be mixed thoroughly immediately before analysis. There must be documentation that specimen mixing by rotary mixer, rocker, automated sampler, or manual inversions is sufficient to ensure reproducibility of CBC results. Some rocking platforms may be adequate to maintain even cellular distribution of previously well mixed specimens, but are incapable of fully mixing a settled specimen. For instruments with automated samplers, the laboratory must ensure that the automated mixing time is sufficient to homogeneously disperse the cells in a settled specimen.

Specimen Collection and Preparation	<ul> <li>Must use fresh blood samples</li> <li>Blood samples must be properly filled to tube manufacturer specifications</li> </ul>
Materials Required But Not Provided	10 fresh normal blood samples

### **Procedure** The following is a suggested method for testing adequate mixing.

- 1. Obtain 10 fresh normal blood samples. Make sure the tubes are filled properly. Do not use under-filled sample tubes.
- 2. Premix the blood samples on a laboratory rocker for at least 5 minutes.
- 3. Immediately load all 10 samples into a cassette(s) starting with position 1 and analyze. Record results for WBC, RBC, HGB and PLT on the log sheet.
- 4. Remove the sample tubes from the cassette(s) and place the tubes upright into a laboratory test-tube rack. Allow the sample tubes to sit UNDISTURBED for 4 hours.
- 5. After 4 hours, place a single tube into position 1 of a cassette and load onto the instrument WITH NO PREMIXING. Allow the analyzer to perform all mixing. Record the results for this analysis on a log sheet.
- 6. Repeat with the remaining 9 sample tubes run each sample tube individually in position 1 of a cassette with no premixing and record the results. Wait until the previous sample's cassette has unloaded before loading the next sample's cassette.

#### Results

### **IMPORTANT**

The responsibility to determine acceptable performance of data lies with the evaluating laboratory.

Compute the absolute difference, percent difference for each sample and the mean absolute difference and mean percent difference.

**Absolute Difference = Premixed result minus Analyzer mixed result** 

### **Percent Difference = (Absolute Difference / Premixed result) x 100**

#### References

Clinical Laboratory Standards Institute. Procedures for the collection of diagnostic venous blood specimens by venipuncture – seventh edition; approved standard GP41. Wayne, PA: CLSI, 2017.

Clinical Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens – sixth edition; approved standard GP42-A6. Wayne, PA: CLSI, 2008.

Clinical Diagnosis and Management by Laboratory Methods, 23rd edition, R. McPherson, M.D, 2016.

Inspection Checklists and Commentary for Laboratory Accreditation, College of American Pathologists, 2017, http://www.cap.org.

#### NOTE

As of 2005, the name of the National Committee for Clinical Laboratory has changed to Clinical and Laboratory Standards Institute (CLSI).

## Verify Adequate Mixing Logsheet

Instrument: \_\_\_\_\_

Date: \_\_\_\_\_

					-	Tec	:h:			
Number	Sample ID	Premixed WBC	Analyzer mixed WBC	Absolute Difference	% Difference		Premixed RBC	Analyzer mixed RBC	Absolute Difference	% Difference
1										
2										
3										
4										
5										
6										
7										
8										
9										
10							-			
Average										
						_				
Number	Sample ID	Premixed HGB	Analyzer mixed HGB	Absolute Difference	% Difference		Premixed PLT	Analyzer mixed PLT	Absolute Difference	% Difference
1										
2										
3										
4										
5										
6										
7										
8										
9										
10										
Average										

Object Name: Hematology Performance Verification Manual (HPVM)	Effective:	09/17/18 (mm/dd/yy)
Doc HENRAHOLOGO Y BERKEMARCE AND HER	Status:	Released

**COMMENTARY** 8

## Manual Differential versus Automated Differential

In the process of correlating an automated differential system with the traditional manual differential, questions or concerns may be raised by laboratory staff, pathologists, or physicians regarding the effectiveness of these methods. The following is a collection of excerpts, from various articles which may answer any questions or concerns.

## I. Manual Diff (General)

Discussion of factors which affect the precision of the manual differential.

## **Comment 1. Dutcher TF<sup>7</sup> states:**

Prior to judging the worth of a routine 100 or 200 cell differential as a screening test, an admission test, or a repeated follow-up test, a number of factors deserve consideration, including the following:

- Technically imperfect smears
- Non-random distribution of cells in smears and non-random search patterns by technologists
- Lack of consensus of definition of cell types particularly bands, polys, monocytes, and atypical lymphocytes
- Intra and inter-observer inconsistencies of cell identification, particularly with bands, polys, monocytes and atypical lymphocytes
- Statistical effect of counting only 100 to 200 cells, including the probability of finding small populations
- Correlation with disease
- Physiology variability of cell classification

It is widely accepted that the distribution of cells in the usual wedge smear is not random. Monocytes tend to congregate at the edge of the smear, lymphocytes in the central area, and neutrophils in the feathered edge. Further, technologists usually begin the differential search in the center of the smear and move toward the edge, but seldom get there because the morphology of the cells becomes less crisp.

Koepke<sup>9</sup> concluded that except for patients with a hematological disease, more than 75 percent of the total variance in neutrophils, lymphocytes, monocytes, and basophils was due to analytical and sampling variability and physiologic changes; the variance was not due to disease. Considering all of the analytic variables involved, it should be apparent that the traditional leukocyte differential is one of the worst procedures offered by the hematology laboratory.

## General diff Comment 2. Rumke CL<sup>4</sup> states:

A blood contains 10 % monocytes. A large number of survey smears are prepared. In each smear n=100 cells are differentiated...the limits of the 95 % prediction interval are 5 % and 16 %. That is, about 95 % of the survey results will lie between 5 % and 16 %.

## General Diff Comment 3. Keopke JA<sup>9</sup> states:

The performance of [routine, 100-cell] differential leukocyte counting by 73 technologists and technicians working in five different laboratories in a large medical center [ Duke University Medical Center, Durham, NC ] was evaluated...

- Good correlation with the reference method CLSI H20<sup>3</sup> : 200 cell differentials performed by each of four referees on four peripheral blood smears, combined for a total 800-cell count] was found for neutrophils [ r = 0.87, comparing 100-cell diff with reference diff ], normal lymphocytes [ r = 0.73 ], and eosinophils [ r= 0.83 ].
- More variability was noted in the estimation of stab neutrophils [ r= 0.67 ], variant lymphocytes [ r= 0.30 ], and monocytes [ r= 0.41].
- In the Clinical Sensitivity analysis of the data, there was a 20 % False normal rate for 'monocytosis' i.e. 20 % of the specimens determined the 800-cell reference differential to exhibit 'monocytosis' (defined as greater tha 12 % monocytes) were classified as 'normal' by the routine 100-cell differential.
- It is widely appreciated that the 100-cell differential count suffers from relatively poor precision as well as the subjectivity inherent in cell classification.

## General Diff Comment 4. Bessman JD<sup>10</sup> states:

The statistical variation... of a 10 % eosinophil count is nearly 100 % of the value... Band counts are similarly unreliable.

Detection of increase bands adds little to diagnosis of infection of inflammation in subjects with neutrophilia and is of uncertain value in subjects with normal white blood cell counts. Increased bands may precede neutrophilia as a sign of infection. However, the quantitation of bands is so variable that several groups have suggested the band elevation adds little to patient management.<sup>3, 12, 14, 15</sup>

## II. Manual Diff (Bands)

Comments on the clinical utility of the 'band count'.

Band comment 1. Clarke G., MD on behalf of the Hematology working group for CCQLM (Canadian Coalition for Quality in Laboratory Medicine). <sup>27</sup> Toward Abandoning The Band - A practice recommendation of the CCQLM:

## **Conclusions and Recommendations:**

Given the widespread availability of automated differential cell counts and the poor analytic and clinical reliability of band neutrophil counts, band neutrophils should no longer be reported as part of the routine manual WBC differential cell count. The laboratory has a responsibility to convey accurate and clinically useful laboratory test results to the clinical community. Band neutrophil counts are not analytically robust or clinically relevant discriminators and should be abandoned.

Where clinical management protocols require assessment for immature neutrophils, qualitative reporting of a left shift based on the presence of circulating myeloid precursors may provide useful information without conveying potentially misleading numeric information.

Band neutrophils should be counted with segmented neutrophils and the two values together reported as a single absolute neutrophil count. Protocols for reporting a neutrophilic left shift as a qualitative result together with red cell morphologic features may be developed by laboratories when a particular clinical need for such information is identified.

# Band comment 2. Novak R, responded to the following questions in CAP Today<sup>15</sup>:

**Q**: How valid is the I/T [ immature/total neutrophil ] ratio?

**A**: Monroe, et al<sup>16</sup> stated that the abnormal values in I/T ratios were found in 51 percent of neonates with confirmed bacterial disease [ not: therefore 49 % had normal I/T ratio ]...but the I/T ratios also were abnormal in 77 percent...of neonates where there was maternal fever but no neonatal disease.<sup>17</sup> Therefore, the I/T ratio does not seem particularly useful. Some authors...have found the I/T ratio to be helpful, but not much more useful than the abnormal total neutrophil count, which, when generated by a hematology analyzer, is much more rapid and accurate. Given the inherent weaknesses that are present in the 100 cell manual differential, I think the absolute neutrophil count, when automated, is a much more precise determination.

## Band Comment 3. Dutcher TF stated in CAP Today<sup>17</sup>:

An acute demand on bone marrow to release neutrophils in response to an infection or inflammation theoretically should cause a rapid depletion of the mature neutrophil reserves, and bands should be released in increasing proportions, producing more immature neutrophils and a 'left shift'. In practice, 100324381e1d31a

the variability inherent in differentiating only 100 or 200 cells can obscure the shift or, contrarily, produce a 'left shift' where none exists.

According to the generally accepted calculations of Rumke, at least 100 cells must be classified to get a reasonably accurate differential. This conclusion is based on the inescapable statistical variability in the counting of cells in a blood smear and is independent of technical factors relative to the quality of the smear or the microscopist's identification of cells.

Rumke's calculations demonstrated that when a perfect technologist examines a perfect smear in a random search pattern of 100 cells and reports 10 percent bands, the patient-care physician can have 95 percent confidence that the true band count is a number somewhere between two percent and 21 percent...Assuming that a 'normal' band count is from zero to 10 percent, many false-positive counts occur with observed bands between approximately 11 percent and 20 percent, whereas many false negatives occur with observed counts from about five percent to 10 percent. The true band count cannot be assumed to be greater than 10 percent until the observed percentage is greater than approximately 25 percent.

It is important for laboratory employees to educate patient-care personnel about the inaccuracies of the differential so that important clinical decisions will not be made on the basis of trivial changes in the differential. Fortunately, patient care physicians must integrate a number of data items before making a clinical decision. They can use the band count if it fits or discard it if it doesn't fit.

## Band comment 4. Savage RA stated in CAP Today<sup>17</sup>:

My initial plan called for this question [ regarding the significance of the 'band count' as it relates to the blood differential count ] to be answered in a 'point-counter-point' format. Dr Dutcher ably presented the 'con-bands' side of the argument. I also asked two authors of recent publications touting the value of band counts... to present their points of view...Neither of the 'pro-band' authors consented to participate in this debate. Thus...I elected to publish Dr. Dutcher's study with a few comments of my own.

Many of the 'pro-band' studies in the literature are seriously flawed by inadequate data or analysis. For example, one study [ Baron MA, Fink HD. Bacteremia in Private Pediatric Practice. Pediatrics, 1980; 66:171-175 ]<sup>18</sup> presented data indicating that the combination of four parameters [ WBC, neutrophil, band, and poly (seg and band) numbers ] predicted...sepsis...in babies three to 24 months old. I recalculated Bayesian statistics using the data in this paper and determined that...the band criterion alone is almost useless, and deletion of the band criterion seems to have little effect on the predicative value of the remaining parameters in detecting bacteremic pediatric patients.

Emergency physician G.P. Young, MD [ CBC or not CBC? That is the question. Ann Emerg Med, 1986; 15:367-371 ]<sup>19</sup> summarized the literature well: "To 'rule in' bacterial infection, most physicians believe that leukocytosis (total WBC>10x10<sup>3</sup>/ cu mm) accompanied by an increase in immature band forms is both specific and sensitive. A review of the WBC literature does not support this conclusion.

## Band Comment 5. Cornbleet PJ and Novak RW stated in CAP today<sup>20</sup>:

Statistical considerations clearly suggest that a 100-cell manual differential yields a value for bands that is not precise. Relatively few cases of infections or inflammatory disorders are seen in which elevation of bands is not accompanied by either neutrophilia or leukocytosis, parameters that are more precisely measured by automated hematology instruments.

Classifications of bands shows poor reproducibility, because it requires establishing discrete borders in what is essentially a continuous process. Criteria for separating band form segmented neutrophils varies among laboratorians... Obviously, the normal range is markedly influenced by the definition.

## Band Comment 6. Wenz B, Gennis P, Canova C, and Burns ER stated<sup>21</sup>:

Leukocyte differentials from 468 emergency room patients were assessed for clinical value. ...in the absence of leukocytes , an elevated band count was instrumental in suggesting admission for only one patient...It is concluded that most leukocyte differentials performed for emergency room patients provide information that is no more clinically significant than that obtained from the medical history, physical examination, and absolute leukocyte count.

## Band Comment 7. Ardron MG, Westengard JC and Dutcher TF stated<sup>22</sup>:

Hospital patients with a documented infection and total leukocyte counts in the reference range were selected for study...Five hundred cell differentials were performed to reduce sampling error. Only one investigator performed all differentials to reduce observer variability.

Applying a threshold value of 15 % for bands, the sensitivity [% true positives ] was 51.1 %, specificity [% true negative ] was 70.0 % and the efficiency [% patients with infection correctly identified by the band count ] was 60.3 %...These are unsatisfactory low statistical values, despite the great care taken to consistently and correctly identified bands. We believe that it is necessary to turn from manual band counts to other more reliable measure for predicting an infections process...Other parameters such as the total leukocyte count, absolute neutrophil count, and temperature have been proven repeatedly to be much better predictors of infection.

## III Manual Diff (Monocytes)

Evaluators speak of the improved accuracy and precision of the monocyte count using an automated differential system.

# Monocyte Comment 1. Gossens W, Van Hove L and Verwilghen RN Stated<sup>23</sup>:

The visual monocyte count has low precision, mainly due to the small number of cells (usually not more than 100-200 leukocytes) identified routinely in this procedure. When 800 leukocytes are counted [ according to the CLSI standard  $H20^3$ ] on a blood sample containing 5-10 % monocytes, a satisfactory precision is reached. Higher Precision should be reached with methods counting a much larger number of events. [ i.e., automated systems ]

It is clear that problems exist in the evaluation of different instruments for counting monocytes. An accurate and reliable reference method is a prerequisite to evaluate this aspect of cell counts. As the visual method is too cumbersome, a different reference method would be useful. Based on the results to this study, it is suggested that the technique using fluorescence-labeled monoclonal antibodies should be regarded as an acceptable alternative. [ This recommendation has been submitted to the International Council for Standardization in Hematology (ICSH) for consideration ]

# Monocyte comment 2. VCS Technology: Monocyte Counting on the Coulter STKS and MAXM.<sup>25</sup>

As part of the ongoing evolution of VCS Technology, recent improvements for the WBC differential of STKS and MAXM instruments have further increased the accuracy of population discrimination. The greatest impact will be seen in the monocyte parameter, primarily because the 2-3 % [ average ] change is more evident in a cell type which typically exhibits a low frequency.

To assess the impact of the changes on the monocyte percentage count, two separate studies [ using flow analysis of MO2/KC56 (CD45/CD14) monoclonal antibodies ] were performed. The improvements for the differential provide excellent agreement on monocyte counts. Monocyte counts were also compared to the manual method [ 800-cell differential, according to CLSI standard H20-A ] ...at each of two sites. Both sites demonstrate closer agreement with manual monocyte counts [ compared to earlier instrument software ]

The normal ranges [for monocytes] that appear in the texts...vary considerably. It is likely that some physicians may question the monocyte values from the automated analyzer. This makes it all the more important to reestablish or verify the reference ranges for normality in current use, and to communicate the meaning of the 'range' to those physicians.

## Monocyte Comment 3. Seaberg R and Cuomo J stated<sup>6</sup>:

Correlation values of the two methods [Coulter STKS and Beckman FACScan flow cytometer] Vs the 800-cell manual count indicated monocyte counts by the STKS and FACScan to be consistently greater than the manual differential results by 2 to 4 percentage points. Possible explanations for this are in the well-known drawbacks of the manual methods: staining quality varies, distributional errors are present on the slide, and observer bias occurs since the morphological classification of some cells is equivocal. Furthermore, the number of cells in a relatively rare population poses statistical problems. These difficulties are reduced in the FACScan and STKS systems. The STKS [VCS] method is fast, efficient, and cost-effective.

GLOSSARY 9

## Glossary

Accuracy of a test is the closeness of the test result to the true value.

**Binomial Envelope** binomial distributions refer to alternative possibilities and their probabilities. The standard error of the binomial distribution is useful in studies where two methods of greatly differing variance are compared, such as manual vs. automated differential.

*Binomial Distribution* is used to generate binomial envelopes applied to the XY scatterplot of data. Usually, results exceeding the envelope are considered inaccurate within a given level of confidence. In our case, we use Dr. Rumke's 95 % confidence limits.

**95 % Confidence Intervals** describes the range within which the statistic is expected to fall in 95 % of the cases were the experiment to be repeated on an identical one expects that 1.0 and zero will fall between the lower and upper limits for the slope and intercept respectively.

*Carryover* is the amount in percent of sample remaining in the system and picked up by the next sample cycled. Low-to-high carryover is the amount of sample with low cell concentrations carried over to sample with high cell concentration, such as diluent to blood. High-to-low carryover is the amount of sample with high cell concentrations carried over to sample with low cell concentrations such as blood to diluent.

*Coefficient of Variation* (CV %) is the standard deviation expressed in terms of the percentage of the mean.

*Correlation Coefficient* (r) is a way to measure the linear relationship between two groups. A value of 1.000 is a one to one relationship.

*Correlation Coefficient Squared* ( $r^2$ ) is a coefficient of determination: how much variability is associated with changes in "X".

**Deming regression** is a linear regression approach invented by the famous QA guru, W. Edwards Deming. Since this calculation assumes that error exists in the data plotted on the X axis, it is a better statistical description of clinical laboratory data than Regular Regression.

*Efficiency* is the percentage of specimens correctly categorized by the automated differential.. (#TP + #TN/Total) x 100

*False Negative* is the percentage of specimens considered normal by the test method (new instrument) and abnormal by the reference method (manual diff).

*False Positive* is the percentage of specimens considered abnormal by the test method (new instrument) and normal by the reference method (manual diff). *Linear regression* is a statistical technique to find the line that best predicts Y from X.

*Mean* is the center of distribution. It is used to describe or summarize a collection of data be defining an average.

*Mean Difference* is the unknown minus the known.

*Measuring range (linearity)* is a set of values of measurands for which the error of a measuring instrument is intended to lie within specified limits.

*Outlier* is a result outside of the established range.

Pass Rate is the percentage of specimens that do not cause a review.

*Precision* is the closeness of results of repeated analysis performed on the same material. Precision is often expressed quantitatively as standard deviations, or coefficients of variation. (also known as repeatability or reproducibility)

*Predictive Value of a Negative Test* is the percentage of specimens that were normal on both the instrument and the manual differential.

*Predictive Value of a Positive Test* is the percentage of specimens that were abnormal on both the instrument and the manual differential.

*Review Rate* is the percentage of specimens that cause a review.

*Patient Population* is the group of individual values about which a generalization is derived.

*Reference interval* is the interval between, and including the lower reference limit to the upper reference limit of the reference population (e.g., 95 % of persons presumed to be healthy or normal.

*Repeatability* is the closeness of results of repeated analysis performed on the same material. Precision is often expressed quantitatively as standard deviations, or coefficients of variation. (also known as precision or reproducibility)

*Reproducibility* is the closeness of results of repeated analysis performed on the same material. Precision is often expressed quantitatively as standard deviations, or coefficients of variation. (also known as repeatability or precision)

*Sensitivity* is the percentage of manual differential abnormals that were also abnormal on the instrument.

*Specificity* is the percentage of manual differential normals that were also normal on the instrument.

*Standard Deviation* (SD) is a measure of the dispersion or variation in the results about the mean.

*True Negative* is the percentage of specimens considered normal (negative) by both the test method (new instrument) and the reference method (manual diff).

*True Positive* is the percentage of specimens considered abnormal (positive) by both the test method (new instrument) and the reference method (manual diff).

*Truth Table* is a technique for evaluation of method sensitivity, specificity, and efficiency.

- 1. CLSI document EP09 User Comparison of Quantitative Clinical Laboratory Methods Using Patient Samples.
- 2. Westgard J and Hunt M. Use and Interpretation of Common Statistical Tests in Method-comparison Studies. Clinical Chemistry, 1973; 19(1): 49-57.
- 3. CLSI document H20 Reference Leukocyte (WBC) Differential Count (Proportional) and Evaluation of Instrument Methods.
- 4. Rumke, CL. The Statistically Expected Variability In Differential Leukocyte Counting. College of American Pathologists CAP Conference, Aspen, CO, 1977; 39-45.
- 5. Galen RS and Gambino SR. Beyond Normality: The Predictive Value And Efficiency Of Medical Diagnosis. John Wiley and Sons, New York, NY, 1975.
- 6. Seaberg C and Cuomo. Assessment Of Monocyte Counts Derived From Automated Instrumentation. Lab Med, 1993; 24 (4): 222-224.
- 7. Dutcher TF: Leukocyte Differentials: Are they worth it? Paper presented at Symposium on Controversies in the Laboratory, Royal Oak, MI, 1984; March (4): 71-87.
- 8. Bull and Korpman: Characterization of the WBC Differential Count. Blood, 1980; 6:411-419.
- 9. Koepke JA et al. A critical evaluation of the manual/visual differential leukocyte counting method. Blood Cells, 1985; 11: 173-186.
- 10. Bessman JD et al. The 'Diff-If'. American Journal of Clinical Pathology, 1988; 90 (2): 763-768.
- 11. Duncan KL and Gottfried EL. Utility Of The Three-Part Differential Count. American Journal of Clinical Pathology, 1987; 88:308-313.
- 12. Mathy KA and Keoke JA. The Clinical Usefulness Of Segment Vs. Stab Neutrophil Criteria For Differential Leukocyte Counts. American Journal of Clinical Pathology, 1974, 958.
- 13. Rainey P. Utility Of The Leukocyte Differential In The Diagnosis Of Bacterial Infection. Pathologists, 1985; 15:11-24.
- 14. Rich EC, Crowson TW and Connelly DP. Effectiveness Of The Differential Leukocyte Count In Case Finding In The Ambulatory Care Setting. Journal of the American Medical Association, 1983; 240:633-636.
- 15. Novak R. Classifying Bands And Segmented Neutrophils. CAP Today, 1995; 1:66.
- 16. Monroe BL et al. Neonatal Blood In Health And Disease. Journal of Pediatrics, 1979; 95:89-98.

- 17. Dutcher TF (author) and Savage RA (ed) Questions and Answers. CAP Today, 1992; 4:65-66.
- 18. Baron MA, Fink HD. Bacteremia in Private Pedriatic Practice. Pedriatics, 1980; 66:171-175.
- 19. Young GP. CBC or not CBC? That is the Question. Annals of Emergency Medicine, 1986; 15:367-371.
- 20. Cornbleet PJ and Novak RW. Classifying Segmented and Band Neutrophils. CAP Today. 1994; 5: 37-41.
- 21. Wenz B, Gennis P, Canova C, and Burns ER. The Clinical Utility Of The Leukocyte Differential In Emergency Medicine. American Journal of Clinical Pathology, 1983: 298-303.
- 22. Ardron MJ, Westgard JC and Dutcher TF. Band Neutrophil Counts Are Necessary For The Diagnosis Of Infection In Patients With Normal Total Leukocyte Counts. American Journal of Clinical Pathology, 1994; 11:646-649.
- 23. Gossens W, VanHove L and Verwilghen RL. Monocyte Counting: Discrepancies In Results Obtained With Different Automated Instruments. Journal of Clinical Pathology, 1991; 44:224-227.
- 24. Westgard, JO. "But is it really out?" Doing the wrong QC wrong. Essay 115. September 2007.
- 25. VCS Technology: Monocyte Counting on Coulter STKS and MAXM. Coulter Corporation. October, 1991.
- 26. Burns, ER. Et al. Evaluation of a Reticulocyte Analysis Method Using a Conventional Automated Blood Cell Differential Instrument. Blood 1982: 10:93a Abstract #358.
- 27. Clarke, G, MD. Toward Abandoning the Band A practice recommendation of the CCQLM. May 2006.
- 28. CLSI Approved Guideline H44, Reticulocyte Analysis Using Supravital Stains and Flow Cytometry.
- 29. COULTER CORPORATION. STKS Analyzer with Reticulocyte Analysis Reference Manual. March 1995.
- 30. Bourke GJ and NcGilvray J. Interpretation and Uses of Medical Statistics. Blackwell Scientific Publications, Osney Mead, Oxford, UK, 1975.
- 31. CLSI document C28 Defining, Establishing and Verifying Reference Intervals in the Clinical Laboratory.
- 32. CLSI document C24 Statistical Quality Control for Quantitative Measurement Procedures: Principles and definitions; 4<sup>th</sup> edition, September 2016, approved guideline.

## CLINICAL LABORATORY STANDARDS INSTITUTE (CLSI) GUIDELINES

#### Standards and Guidelines:

- **EP 09 A2** Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline Second Edition (2011)
- **EP 10 A3** Preliminary Evaluation of Quantitative Clinical Laboratory Measurement Procedures; Approved Guideline – Third Edition (2014)
- **EP 15 A3** User Verification of Precision and Estimation of Bias; Approved Guideline Third Edition (2014)
- **GP39** Tubes and Additives for Venous Blood Specimen Collection Fifth Edition; Approved Standard Sixth Edition (2010)
- GP41 Collection of Diagnostic Venous Blood Specimens Seventh Edition (2017)
- **GP42** Procedure and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard – Sixth Edition (2008)
- **GP44** Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline Forth Edition (2010)
- **H20 A2** Reference Leukocyte Differential Count (proportional) and Evaluation of Instrumental Methods; Approved Standard Second Edition (2007)
- **H26** Validation, Verification, and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition (2010)
- **H44 A2** Methods for Reticulocyte Counting (Automated Blood Cell Counters, Flow Cytometry and Supravital Dyes); Approved Guideline (2004)
- H56 A Body Fluid Analysis for Cellular Composition; Approved Guideline (2006)

Order by Phone	Toll-free	(877) 447-1888
Order by Fax	(610) 688-070	00
Website	www.clsi.org	
E-mail	customerservi	ice@clsi.org

## WEBSITES FOR LABORATORY REFERENCE

Beckman Coulter Inc. www.beckmancoulter.com

Clinical Laboratory Standards Institute (CLSI) www.clsi.org

College of American Pathologist (CAP) www.cap.org

Joint Commission (JCAHO) www.jointcommission.org

Centers for Medicare and Medicaid Services (CMS) www.cms.hhs.gov

Clinical Laboratory Improvements Amendments www.cms.hhs.gov/clia

- Proficiency Testing Providers www.cms.hhs.gov/CLIA/14\_Proficiency\_Testing\_Providers.asp
- Media Lab, Inc. (Online Compliance and Continuing Education for Clinical Laboratories www.medialabinc.net

Westgard QC (Tools, Technologies and Training for Healthcare Laboratories) www.westgard.com

Lab Tests Online www.labtestsonline.org

Object Name: Hematology Performance Verification Manual (HPVM) PERFORMANCE VERIFE ACTION MANUAL (Mmm/dd/yy) Doc Title: Hematology Performance Verification Manual (HPVM) Status: Released

P/N: 4277076AC