



# Maintenance Checklist

CellaVision DM1200

Serial # \_\_\_\_\_

System ID \_\_\_\_\_

Month \_\_\_\_\_ Year \_\_\_\_\_ 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 24 25 26 27 28 29 30 31

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
<b>DAILY</b>																																
QC: Cell Location - (≥ 97%)																																
Perform BF Cell Location - as applicable																																
<b>Tech Initials</b>																																
<b>WEEKLY</b>																																
Clean Objectives and LED Table																																
Clean Hood																																
Clean Bottom Tray																																
Power off/on System Computer																																
Check Immersion Oil Level																																
Delete unsigned and failed orders																																
<b>Tech Initials</b>																																
<b>REMEDIAL</b>																																
Change Immersion Oil Pack																																
Remove Stuck Slide from the Gripper																																
<b>DATABASE PERFORMANCE</b>																																
Check the Database Size																																
Activate Autodelete/Archiving, if Necessary																																
<b>Tech Initials</b>																																