



DxH 520 Instrument Maintenance Checklist

Month _____ Year _____

Serial # _____

System ID _____

DAILY	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Perform Shutdown																																
Perform Daily Checks - Verify that all results are within limits																																
Clean the Instrument																																
Process Controls																																
Monthly or every 1,000 cycles																																
Perform Bleach Cycle																																
Monthly																																
Clean the WBC Bath Filter																																
As Needed																																
Clean the Tube Holder and Tube-Holder Housing																																
Clean the Baths																																
Clean the Bar-code Scanner																																
Drain the Baths																																
Backflush the Apertures																																
Yearly																																
Lubricate the Pistons																																
Yearly or every 18,000 Cycles																																
Replace the Rinsing Head O-Ring																																
Replacement As Needed																																
Replace the Aspiration Probe																																
Remove or Replace the Tube Holder																																
Replace the Bar Code Scanner																																
Tech Initials																																