

OPERATOR TIPS

UniCel® DxC 600/800 SYNCHRON® Clinical System
 STORE **UN-OPENED** CALIBRATORS & CONTROLS AS FOLLOWS:

CALIBRATORS / CONTROLS

ROOM TEMPERATURE	REFRIGERATOR 2°C to 8°C	FREEZER -15°C to -20°C
<i>CALIBRATOR</i>	<i>CALIBRATOR</i>	<i>CALIBRATOR</i>
FE/IBCT Calibrator	Ammonia Calibrator	APO Calibrator
	CAL 5 Plus	Bilirubin Calibrator
	CRP Calibrator	Enzyme Validator
	VANC Calibrator	Lipid Calibrator
	DAT Calibrators	M-TP Calibrator
	DRUG Cal 1	Multi Calibrator
	DRUG Cal 2	Protein Calibrator
	DRUG Cal 3 Plus	
	ETOH Calibrator	
	HbA1c/HbA1c2 Calibrator	
	HDL Calibrator	
	LDLD Calibrator	
	Lipase Calibrator	
	Lithium Calibrator	
	LX/DxC Cal 1	
	LX/DxC Aqua Cal 1, 2, 3	
	LX/DxC MA Calibrator	
	LX/DxC PAB Calibrator	
	Salicylate Calibrator	
	T4 Calibrator	
	T-Uptake Calibrator	
	<i>CONTROLS</i>	<i>CONTROLS</i>
	ALC/AMM Controls	All Vigil Controls
	DAT Controls	NYSPATH Controls
		SYNCHRON Controls
		TRIAD Controls
		Ultimate Bilirubin Controls