



UniCel DxH Series Workcell Cleaning Checklist

Month _____ Year _____

Serial # _____

System ID _____

DAILY	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	
Perform Shutdown - Workcell																																
Perform Daily Checks - Verify that all results are within limits - Workcell																																
Process Controls - DxH 800																																
Cleaning Stainer Baths and Tray (Daily or as Needed) - DxH SMS																																
Weekly or As Needed																																
Cleaning the Baskets - DxH SMS																																
Clean Stainer Fill Probes, Drain Probes and Level - DxH SMS																																
Flush Stainer Module (Automatic) - DxH SMS																																
Flushing Reagents Lines (Manual) - DxH SMS																																
Extensive Cleaning of Fill Probes and Drain Probes - DxH SMS																																
As Needed																																
Clean the BSV Externally - DxH 800																																
Clean (Bleach the Aperture) - DxH 800																																
Clean the Pneumatic Supply Module Fan Filter - DxH 800																																
Clean the (AMTC) Module - DxH 800																																
Cleaning the STM - Workcell																																
Cleaning the Aspiration Probe - Workcell																																
Cleaning the STM Optical Sensors - Workcell																																
Cleaning the Vacuum Trap - Workcell																																
Cleaning the STM Bar Code Scanner - Workcell																																
Loading Slides - DxH SMS																																
Cleaning the Dispense Probe - DxH SMS																																
Emptying the Broken Slides Bin - DxH SMS																																
Cleaning the Slide Chute - DxH SMS																																
Tech Initials																																