



UniCel DxH Slidemaker Stainer Cleaning Checklist

Month _____ Year _____

Serial # _____

System ID _____

DAILY	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Perform Shutdown																															
Cleaning Stainer Baths and Tray (Daily or as Needed)																															
Weekly or As Needed																															
Cleaning the Baskets																															
Clean Stainer Fill Probes, Drain Probes and Level																															
Flush Stainer Module (Automatic)																															
Flushing Reagents Lines (Manual)																															
Extensive Cleaning of Fill Probes and Drain Probes																															
As Needed																															
Cleaning the STM																															
Cleaning the Aspiration Probe																															
Cleaning the STM Optical Sensors																															
Cleaning the Vacuum Trap																															
Cleaning the STM Bar Code Scanner																															
Loading Slides																															
Cleaning the Dispense Probe																															
Emptying the Broken Slides Bin																															
Cleaning the Slide Chute																															
Tech Initials																															